

PLAS-ILABS™

Clearly Your Best Choice™



CONCEPTION ■ INTERPRETATION ■ PRODUCTION ■ SATISFACTION

PLAS ■ LABS™



*Quality research
like quality music...*

*Should begin with
fine equipment.*

*When you are
ready to start your
research in grand
style, contact the
creative people
at Plas ■ Labs, Inc.*

*We are approaching
our fifth decade
of manufacturing
specialized research
equipment. Chambers
are one of the keys
to our success.
Can we get in tune
with yours?*

PLAS ■ LABS™

401 E. NORTH ST.
LANSING, MICHIGAN
48906 U.S.A.

Innovated Quality Since 1967

www.plas-labs.com ■ E-mail: sales@plas-labs.com ■ 1.517.372.7177 ■ 1.800.866.7527 ■ Fax: 1.517.372.2857

OUR GLOBAL DISTRIBUTION ARRANGEMENTS

PLAS ■ LABS Inc.

Our Global Distribution Network.

Company established in 1967.

Plas-Labs has a global distribution network of over 110 distributors. Our distributors have a proven record of good service, applications expertise, and knowledge of our product line. Many are scientists and engineers with strong technical backgrounds and “hands-on” experience.

Most of our global distributors are family owned businesses. Good service and customer satisfaction are critical to their ongoing business success.

For a complete listing of our distributors (country by country) please go to the Plas-Labs web site at www.plas-labs.com

Plas-Labs has a global pricing policy maintained on a year long basis. We normally change prices once a year. All quotations are honored for sixty (60) days.

We are proud of our products, their service and durability. Our product warranty is twenty-four (24) months, but not including gloves and consumables.

Have a need or an idea for a new product?

Just let us know.

We like to do custom work !

Phone: 1-517-372-7177

Fax: 1-517-372-2857

“e” sales@plas-labs.com



“CLEARLY YOUR BEST CHOICE”™

WHAT'S INSIDE AND WHAT'S NEW !! • INDEX

Automatic Controlled Atmosphere (Anaerobic) Chamber (<i>Touch Screen Controlled</i>)	4	NEW!
Automatic Controlled Atmosphere (Anaerobic) Chamber (<i>"PLC" Controlled</i>)	5	NEW!
Controlled Atmosphere (Anaerobic) Chamber	6-7	
Multiple Operator Atmosphere (Anaerobic) Chamber (<i>All Manually Controlled</i>)		
Closed Loop Containment Isolator Glove Box	8-9	
Nitrogen Dry Box	10-11	
Multiple Operator Nitrogen Dry Box (Manually Controlled)		
Basic Glove Box, Single Operator	12-13	
Multiple Operator, Basic Glove Box		
Compact Glove Box (#830-SERIES)	14-15	
Compact Glove Boxes with HEPA Filtration	16	
Compact Glove Boxes with Access Doors & Automatic Locks	17	NEW!
Germicidal Work Station (Still Air Enclosure)	18	
La Petite Glove Box Re-engineered Style!	19	NEW!
PCR (UV) Chamber	20-21	NEW!
HEPA Filtered, PCR Chamber		
Low Velocity Fume Hoods	22-23	NEW!
Isolation & Containment Chambers	24	
Isolation & Sterility Chambers	25	NEW!
Biotoxin & Decontamination Chambers	26	NEW!
Custom Unit Capabilities (Containment & Drug Processing)	27	
Desiccator Cabinets (Clear and Light Sensitive Units)	28-29	
CO ₂ Lab Chambers (Heavy Duty)	30	
Stainless Steel Desiccators	31	
Stainless Steel Glove Box (Integrity) Re-engineered Product!	32	
Vacuum Chambers (Multi Purpose)	33	
ACCESSORIES FOR ALL CHAMBERS AND GLOVE BOXES		
Major Accessories (Devices)	34	NEW!
Digital Oxygen Indicator		
"Counter Spi" Static Eliminator		
Humidity Control System		
Major Accessories (Devices)	35	NEW!
Controlled Purge System		
Bacti-Microcinerator		
"Bright Light" System		

WHAT'S INSIDE AND WHAT'S NEW !! • INDEX

HEPA Filtration Components & Variations	36	
HEPA Filters		
Carbon Filters		
Magnehelic Gauges		
HEPA Filtration Components & Variations (Continued)	37	
HEPA Single System		
HEPA Double System		
HEPA Open Exhaust System		
Suspended Shelf (Storage) Systems	38	NEW!
#830-SERIES Units		
Glove Boxes & Isolators		
Entry Port (Plug) Variations	39	
Bare Hand Entry		
Atmosphere Purification & Sterilization	40-41	NEW!
Carbon & General Filter Types		
“VHP” Sterilization Set Ups		
Catalyst Heater & Palladium Canister		
Drying Train Set Up		
“ONEG” Strict Oxygen Reducing System		
General Purpose Vacuum Pump		
Stainless Steel Support Carts (Trolleys)	42	
“Knocked Down Style”		
Microscope & Inspection Lens	43	NEW!
Glove Box & Isolator Gloves	44-45	
Glove Specifications		
Valve & Gauge Variations	46	NEW!
Banjo™ Fittings		
Magnehelic Gauge Set Ups		
Miscellaneous Accessories	47	NEW!
Blanking Plates		
Computer Interfaces		
Automatic Door Inter Locking System		
Miscellaneous Accessories (Continued)	48	NEW!
Germicidal Dip (Entry) Tank		
“T” Connecting Transfer Chamber		
Oxygen Indication (Hand Held)		

Do you have technical questions and or suggestions? Please contact us at....

Phone: 1-517-372-7177 Fax: 1-517-372-2857 “e” sales@plas-labs.com

CHAMBERS AND GLOVE BOXES

TOUCH SCREEN, AUTOMATICALLY CONTROLLED ATMOSPHERE ISOLATOR/GLOVE BOX

"NEW PRODUCT"



The NEW **#890-THC** accurately maintains temperature and humidity (Rh) inside a roomy single operator glove box. We also offer multiple operator units.

The unit is operated by merely touching the color control screen for the functions desired. It includes your private password protected administrator controls. Touch screen viewing of the data logging function is included as well as serial port connection for data logging.

Automatic purging cycles can be programmed.
The humidity range control accuracy is $\pm 2\%$ Rh.

Temperature can be regulated in either "F" or "C" degrees and the Peltier style thermoelectric assembly is accurate to ± 2 degrees C. Temperature range is from 16 to 44 degrees C.

Features included are:

- Password protected, configurable "PID" controlled.
- Inert gas hookup capability.
- White ambidextrous Hypalon™ gloves.
- Platinum RTD with sealing gland.
- Four (4) pre-set conditions.
- 24 hour one touch, on screen data logging.
- Room temperature and Rh read outs.
- Thermoelectric assembly (Peltier device).
- Nebulizer for the Rh system.
- The top section is one piece clear acrylic plastic for easy viewing. It features "Easy Clean" corners.
- Four (4) key cock valves for purging. Two (2) on the main chamber and two (2) on the transfer chamber.
- 8" (203 mm) Color panel touch screen. Date and time display.
- Humidity (Rh) transmitter with 8 feet of sensor cable.
- Standard 24 month warranty (Not gloves or consumables).
- Serial port for external computer data logging & remote access.
- Transfer chamber is 11" (280 mm) I.D. x 12" (305 mm) long.
- "Bright Light" system (40,000 hour lamp guarantee).

Optional Drying Train package. Catalog #890-THC/DT

The Drying Train accessory includes two (2) canisters filled with Molecular Sieve™ and circulation (vacuum) pump. This accessory can be used if a house supply of inert gas is not available. Periodic regeneration of the Molecular Sieve™ is recommended for optimal drying performance. Refer to Accessories Section.

NOTE: If drying train feature is preferred, order as:
#890-THC/DT (Domestic) or **#890-THC/DT/EXP** (Export).



#890-THC & #890-THC/EXP

Product Specifications

Inside Dimensions	41" w x 28" d x 26" h 1,040mm x 712mm x 660mm h
*Outside Dimensions	55" w x 35" d x 38" h 1,400mm x 890mm x 965mm h
Approximate Volume	17.3 Cubic Feet 489 Liters
Electrical Requirements	N. American 115-120 Volts, 60 Hz. 10 Amps Export Model 220-240 Volts, 50 Hz. 5 Amps.
Maximum Working Pressure	6.0" of W.C. (1.49kPa)
Approximate Shipping (crated) Weight	250 Pounds 112.5 Kilos

*Includes Transfer Chamber



Close up of computer interface screen

CHAMBERS AND GLOVE BOXES

AUTOMATIC “PLC” CONTROLLED ATMOSPHERE ISOLATOR/GLOVE BOX (ANAEROBIC CHAMBER)

This chamber has been designed for automatic atmosphere control when working with Oxygen sensitive materials, product containment needs, and/or general isolation control. The chamber features an internal “PLC” control system with multifunctional operator interface.

This glove box is perfect for working in anaerobic conditions.

Features included are:

- Internal “PLC” control system.
- Three (3) drying train canisters filled with Molecular Sieve™ and quick dis-connects.
- White ambidextrous Hypalon™ gloves.
- Two (2) 1/4” hose barb (tubing) gas connections.
- “Bright Light” system (40,000 hour lamp guarantee).
- Programmable Catalyst Heater unit with two (2) Palladium (Pd) filled canisters.
- Multiple electric outlet (socket) strip. UL, CSA, and CE approved.
- Two (2) ground key cock valves for purging the main chamber.
- Standard 24 month warranty (Not gloves or consumables).

The operator can program these features:

- The number of purging cycles.
- Gas selection (inert or anaerobic blend).
- Drying train reaction time frame.
- Pressure “hold” function.
- Anaerobic gas reaction time frame.
- Selection of chamber to be purged (Main or transfer).

“PLC” Control Board close up



The operator interface continually displays the “PLC” functions and working pressures within the main chamber or transfer chamber. It can also support an additional 4-20 mAmp analog input. These inputs are used to create your own unique atmosphere chamber. Examples are percentage of O₂ analyzer, humidity information, etc.



SIMPLICITY-888 & SIMPLICITY-888/EXP With Catalyst Heater Unit and Palladium Canisters

Product Specifications

Inside Dimensions	41" w x 28" d x 26" h 1,040 mm x 712 mm x 660 mm h
*Outside Dimensions	55" w x 35" d x 38" h 1,400 mm x 890 mm x 965 mm h
Approximate Volume	17.3 Cubic Feet 489 Liters
Electrical Requirements	N. American 115-120 Volts, 60 Hz. 10 Amps. Export Model 220-240 Volts, 50 Hz. 5 Amps.
Maximum Working Pressure	6.0" of H ₂ O 1.49kPa
Approximate Shipping (crated) weight	250 Pounds 112.5 Kilos

*Includes transfer chamber.



#800- HEATER w/#800- PC

The programmable Catalyst Heater has two basic functions. It reduces trace amounts of Oxygen that enter the main chamber and it maintains incubation temperatures from ambient up to 42 degrees C.

*Refer to Accessories Section for more details about this device.

CHAMBERS AND GLOVE BOXES

CONTROLLED ATMOSPHERE (ANAEROBIC CHAMBER) GLOVE BOX

The #855-AC (Anaerobic Chamber) has been designed for manually controlling your atmosphere when working with Oxygen sensitive materials. It is perfect for working in situations requiring low levels of Oxygen during critical isolation of the operator or your research materials. Any inert gas may be used.

Its main advantage is minimizing the risk of Oxygen entry into the main working chamber. If trace amounts of Oxygen happen to enter the main chamber, the Catalyst Heater reduces it to water vapor. The drying train can then be activated to remove any excess water vapor.

The #855-AC is offered as a complete system ready to use. There are no other components required. All you need to add are your work samples and gas of choice.

This glove box has a capacity of 348 culture plates when using our suspended shelf accessory. Please refer to the Accessories Section (#800-SHELF-1).

Some typical applications include:

Microbiology, Biochemistry, Plasma environment work, animal science studies, and electronic sub-assembly work.

The Basic Chemical Reaction.

The internal atmosphere is drawn off the floor of the glove box and is pushed up through the heating elements and the warm Palladium (Pd) pellets. The Hydrogen in the special gas mixture then reacts with trace amounts of Oxygen to create water vapor. The excess water vapor is then absorbed by the Molecular Sieve™ in the three drying train canisters. The canisters are re-chargeable.

The Catalyst Heater Function.

The catalyst heater has two main functions. It reduces trace amounts of Oxygen and it maintains the correct incubation temperature level for cultivating anaerobes (35 to 39 degrees C).

Did you Know?
Molecular Sieve™
can absorb 20-25% of
its own weight in
water (H₂O).



#855-AC with Catalyst Heater Unit
#855-AC/EXP (Export Model)

Refer to #800-ONEG and dimensions on page #7.

Features included are:

- Standard 24 month warranty (Not gloves or consumables).
- Formed one-piece clear plastic top section with “Easy Clean” corners.
- Matched die molded white thermoset bottom with “Easy Clean” corners.
- Two vacuum diaphragm pumps, one each for the drying train and the transfer chamber (purging).
- Drying train includes three (3) clear plastic canisters filled with Molecular Sieve™.
- All controls are illuminated.
- “Bright Light” system (40,000 hour lamp guarantee).
- White ambidextrous Hypalon™ gloves.
- All clamps are adjustable to compensate for wear.
- Adjustable vacuum gauge on transfer chamber.
- Four (4) ground key cock valves for purging.
- Transfer chamber is 12” (305 mm) long x 11” (280 mm) I.D.
- Electrical outlet (socket) strip (UL, CSA, & CE).
- Self sealing quick disconnects allow changing of the drying train without disturbing the internal atmosphere.

Electrical requirements.

North American 115-120 Volts, 60 Hz. 10 Amps.
Export Models 220-240 Volts, 50 Hz. 5 Amps.

CHAMBERS AND GLOVE BOXES

CONTROLLED ATMOSPHERE (ANAEROBIC) CHAMBER FOR MULTIPLE OPERATORS

Our #855-ACB contains all the features of our single operator unit, but it has more capacity. It is ideal for two or more operators working together up to a maximum of four (4) people. Two people each side.

NOTE: The maximum working pressure is 6" of WC.

The Drying Train Function.

The drying train includes its own vacuum pressure pump and six (6) polycarbonate canisters filled with Molecular Sieve™. The canisters are attached in series on the top and back of the unit. All canisters are connected using self sealing quick disconnects so as not to disturb the inner atmosphere when recharging them.

The Molecular Sieve™ acts as an atmosphere purifying medium. Most impurities produced by research will adhere to the granules. It also acts as a drying agent as moisture is produced. As moisture is absorbed, the Molecular Sieve™ will turn color to a dark brown.

The Molecular Sieve™ can be rejuvenated by heating in a hot oven for three (3) hours at 450° F.

NOTE: You must remove the Molecular Sieve™ from the plastic canisters before recharging (baking).

What's included with the glove box:

- Two (2) pairs of white ambidextrous Hypalon™ gloves.
- Two (2) vacuum pressure pumps. One each for the drying train and the transfer chamber purging cycle.
- Two (2) Catalyst Heaters and four (4) Palladium canisters.
- Manual controls for purging, lights, and drying train operation.
- A "Bright Light" illumination system with a 40,000 hour lamp guarantee.
- Hospital grade multiple electric outlet (socket) strip.
- Optically clear plastic top and bottom sections with "Easy Clean" corners.
- Transfer chamber size 12" long X 11" I.D. (305 mm X 280 mm I.D.).
- Two year (24 month) standard warranty (Not gloves or consumables).

Electrical requirements:

North American 115-120 Volts, 60 Hz. 15 Amps.
Export models 220-240 Volts, 50 Hz. 8 Amps.




#855-ACB for Multiple Operators

#855-ACB/EXP for Export Use

Product Data	#855-AC	#855-ACB
I.D.	41"w x 28"d x 26"h 1,040 mm x 710 mm x 660 mm h	60"w x 38"d x 31"h 1,520 mm x 960 mm x 790 mm h
O.D. Includes Transfer Chamber	55"w x 35"d x 38"h 1,400 mm x 890 mm x 970 mm h	76"w x 47"d x 42"h 1,930 mm x 1,190 mm x 1,070 mm h
Approximate Volume	17.3 Cubic Ft. 489 Liters	40.9 Cubic Ft. 1,157 Liters
Approximate Shipping Wt. (Crated)	250 Pounds 113 Kilos	300 Pounds 136 Kilos

#800-ONEG Oxygen reducing accessory w/pump not included

This accessory further reduces the trace amounts of Oxygen remaining in the glove box after initial purging. Part #1 is the



Oxygen trap. Part #2 is the Oxygen indicator (sensor). The accessory is attached only after the initial purge cycles have been completed. It will reduce trace amounts of Oxygen down to about ten (10) ppb.

Refer also to Page #41.

CHAMBERS AND GLOVE BOXES

CLOSED LOOP CONTAINMENT ISOLATOR FOR MULTIPLE OPERATORS

The **#870-CLC** “Closed Loop” Isolator is a completely automated negative pressure containment glove box.

The system is engineered and programmed to maintain a slight negative working pressure to ensure containment of your specimens (-0.5” of H₂O). The internal atmosphere is continually filtered by two (2) HEPA filters. Filters are rated at 99.9997% efficiency and can be changed over from within the isolator for safety. During change over, contaminated filters do not come into direct contact with the operator.

These isolators may be coupled together in train for multiple station work processes. They can be set up for up to four (4) sets of glove ports (operators). Two operators on both the front and back walls.



#870-CLC with Stainless Steel Support Cart
#870-CLC/EXP Export Model with Trolley

Basic features of the CLC Isolator.

- Effective containment of potentially hazardous toxins, organisms, and other materials etc.
- A completely sealed containment unit with automatic negative pressure control.
- Large, easy to read control panel with digital controllers and bright LED status lights.
- Easy (and safe) to change, inlet and outlet HEPA filters. Filters can be incinerated after use.
- Easy to clean highly radiused corners.
- Stainless steel support cart (trolley) with gas tank support tracks for safety.
- Clear plastic top and bottom section. Bottom is formed of high impact resistant plastic.
- “Bright Light” illumination system with a 40,000 hour lamp guarantee.
- Two (2) pairs of white ambidextrous Hypalon™ gloves.
- Manual house exhaust valve for emergency venting of the isolator.
- Electrical outlet (socket) strip (UL, CSA, & CE).
- Four (4) ground key cock valves for purging.
- Standard twenty-four (24) month warranty (Not gloves or consumables).

Product Details	#870-CLC & #870-CLC/EXP
I.D.	56”W x 35”D x 27” high 1,423 mm x 890 mm x 690 mm h
O.D. (Includes Transfer Chamber)	82”W x 47”D x 53” h 2,080 mm x 1,190 mm x 1,350 mm h
Approximate Volume of Isolator.	40 Cubic Feet 1,157 Liters
Approximate shipping weight of isolator & Cart. Two (2) Crates.	Isolator @ 300 Lbs. 136 Kilos. Cart (Trolley) @ 140 Lbs. 64 Kilos

Electrical requirements.

North American 115-120 Volts, 60 Hz. 10 Amps.
Export Models 220-240 Volts, 50 Hz. 5 Amps.

NOTE: Activated Carbon (charcoal) filters are also available for this unit. Four (4) different compositions can be used. Refer to Catalog Accessories Section and contact us for details.

CHAMBERS AND GLOVE BOXES



#870-CLC & #870-CLC/EXP (Close up of top shroud)

Working
with fine powders and
chemicals?

The #800-AS/SPI Ionizer emits both negative and positive ions. The non air assisted static eliminator does not cause any air turbulence because there is no internal fan (blower). No particle contaminant's are introduced into the isolator. It is perfect for working with fine powders and chemicals.


Refer also to page #34.



#800-AS/SPI (Anti-Static unit)
#800-AS/SPI/EXPORT Model

Electrical requirements.

North American 115-120 Volts, 60 Hz. 1 Amp.
Export Models 220-240 Volts, 50 Hz. .5 Amps.



#870-CLC Atmosphere Flow Schematic (Front view)

How the isolator works:

If the internal negative pressure rises above -0.55" of W.C, the programmable logic controller automatically opens the atmosphere supply valve. Incoming (make-up) gas is automatically filtered through the inlet HEPA filter.

Once the internal atmosphere reaches your preferred set point, the automatic valve closes. The isolator also includes a manual "ball valve" for connection to your house exhaust/filter system.

CHAMBERS AND GLOVE BOXES

NITROGEN DRY BOX

Our **#850-NB** "Nitrogen Dry Box" is a completely enclosed environmental chamber ideal for working in your preferred inert gas atmosphere.

Any inert gas such as Nitrogen, Argon, Helium and Plasma are easily maintained.

The unit is compact, lightweight, and portable. The design makes it easily moveable from lab to lab.

Everything you need is included with the glove box. All you need to add is your gas of choice and your preferred research materials.

Basic Features that are included:

- Standard 24 month warranty (Not gloves or consumables).
- Formed one piece clear acrylic top section with rounded "Easy Clean" corners.
- Matched die molded white thermoset bottom with rounded "Easy Clean" corners.
- Two (2) vacuum diaphragm pumps. One each for the drying train and transfer chamber.
- The drying train includes three (3) clear acrylic canisters filled with Molecular Sieve™.
- Transfer chamber is 12" (305 mm) long x 11" (280 mm) I.D.
- All controls are illuminated.
- "Bright Light" system (40,000 hour lamp guarantee).
- White ambidextrous Hypalon™ gloves.
- All clamps are adjustable to compensate for wear.
- The front sloped viewing panel minimizes glare.
- Vacuum gauge on transfer chamber.
- Self sealing quick disconnects enable changing of the drying train canisters without disturbing the internal atmosphere.
- Four (4) ground key cock valves for purging.
- Transfer chamber leveling tray for handling liquids.
- Hospital grade multiple electric outlet (socket) strip. UL, CSA, and CE approved.

Did you Know?
Molecular Sieve™ can be recharged in a hot oven at 450°F for 3 to 4 hours. You must remove the sieve from the plastic canisters before baking however.



#850-NB with optional shelf accessory*
#850-NB/EXP for Export use

***800-SHELF-I**, shown as an optional accessory.
Culture dish capacity from 348 up to 432 plates.
Refer to Accessories Section for more details.

Product Specifications

Product Data	#850-NB	#850-NBB
I.D.	41"w x 28"d x 26"h 1,040 mm x 710 mm x 660 mm h	60"w x 38"d x 31"h 1,520 mm x 960 mm x 790 mm h
O.D. Includes Transfer Chamber	55"w x 35"d x 38"h 1,400 mm x 890 mm x 970 mm h	76"w x 47"d x 42"h 1,930 mm x 1,190 mm x 1,070 mm h
Approximate Volume	17.3 Cubic Ft. 489 Liters	40.9 Cubic Ft. 1,157 Liters
Approximate Shipping Wt. (Crated)	250 Pounds 112.5 Kilos	300 Pounds 135 Kilos

Electrical requirements:

North American 115-120 Volts, 60 Hz. 10 Amps.
Export models 220-240 Volts, 50 Hz. 5 Amps.

NOTE: If this unit is to be used for Anaerobic work, you must also purchase our **#800-HEATER** and **#800-PC** (Palladium) accessories for optimum performance. Refer to the Accessories Section for details.

CHAMBERS AND GLOVE BOXES

NITROGEN DRY BOX FOR MULTIPLE OPERATOR USE

The **#850-NBB** is our Nitrogen Dry Box for two (2) operators. It can be set up for up to four (4) operators. It includes all of the same features as our **#850-NB**, but it has six (6) drying train canisters and two (2) vacuum diaphragm pumps.

NOTE: The unit as shown includes optional accessories **#800-SHELF-II** and **#CART-GBB**, and **#800-CPS/O2**.

Working with fine powders and chemicals?


The **#800 AS/SPI** "Counter SPI" ionizer emits alternating positive and negative Ions that eliminate static electricity within the glove box. It is non-air assisted (no fan), so there is no atmosphere turbulence inside the chamber. Additionally, no particle contaminants are introduced because of the ionizer.

This device is perfect for working with fine powders and chemicals.



#800-AS/SPI Anti-Static unit

Electrical requirements for **#800-AS/SPI**.
North American 115-120 Volts, 60 Hz. 1Amp.
Export models 220-240 Volts, 50 Hz. .5 Amp.



#850-NBB shown with optional accessories
#850-NBB/EXP for Export Use

Optional accessories shown above include:

#800-CART/GBB Refer to page #42

#800-SHELF-II Refer to page #38

#800-CPS/O2 Refer to page #35
(Oxygen Controlled Purge Unit)

For Complete Product Specifications
refer to page #10 (**#850-NB**)

NOTE: If the **#850-NBB** glove box is used for anaerobic work, you must purchase additional accessories **#800-HEATER** and **#800-PC**. For the large (NBB) glove box, we recommend **two (2) #800-HEATERS** and **four (4) #800-PC canisters** for the best efficiency.

CHAMBERS AND GLOVE BOXES

BASIC GLOVE BOXES



Our catalog #818-GB style glove boxes are designed to be just what they are named, “Basic.” This is a simple enclosure that can be used for general isolation or containment work. They are very useful for isolating sensitive research materials from a hostile exterior environment.

These units are engineered to fit general laboratory research situations. They can be easily modified for specific research projects.

These are “starter” glove boxes upon which many accessories can be attached allowing for specific research projects. For unique research applications, please refer to Page #11 and our “Accessories Section” of this catalog.

These units are compact, portable, light weight, and can easily be moved from lab to lab.



#818-GB & #818-GB/EXP Basic Glove Box

Product Specifications

Product Data	#818-GB #818-GB/EXP	#818-GBB #818-GBB/EXP
I.D.	41”w x 28”d x 26”h 1,040 mm x 710 mm x 660 mm h	60”w x 38”d x 31”h 1,520 mm x 960 mm x 790 mm h
O.D. Includes Transfer Chamber	55”w x 35”d x 27”h 1,400 mm x 890 mm x 685 mm h	76”w x 47”d x 32”h 1,930 mm x 1,190 mm x 812 mm h
Approximate Volume	17.3 Cubic Ft. 489 Liters	40.9 Cubic Ft. 1,157 Liters
Approximate Shipping Wt. (Crated)	250 Pounds 113 Kilos	300 Pounds 135 Kilos

Basic Features that are included:

- Standard 24 month warranty (Not gloves or consumables).
- Formed one piece clear acrylic top section with rounded “Easy Clean” corners.
- Matched die molded white thermoset bottom with rounded “Easy Clean” corners.
- White ambidextrous Hypalon™ gloves.
- All clamps are adjustable to compensate for wear.
- The front sloped viewing panel minimizes glare.
- Vacuum gauge on transfer chamber.
- Four (4) ground key cock valves for purging.
- Transfer chamber leveling tray for handling liquids.
- The transfer chamber is 12” (305 mm) long x 11” (280 mm) I.D.
- Hospital grade multiple electric outlet (socket) strip. UL, CSA, & CE approved (Domestic 12 Amp rating).

Did you Know?

The acceptable leak rate (on volume) for glove boxes is .05% per hour. This is an established standard by The American Glove Box Society.

Electrical requirements:

North American 115-120 Volts, 60 Hz. 12 Amps.
Export models 220-240 Volts, 50 Hz. Amps will vary as to the end users country.

NOTE: Clear PVC or PETG units are also available for “VHP” sterilization needs. Two (2) “Banjo™” type attachment fittings are included and mounted on the glove box. This is a special quote item. Refer to pages #40 and #46.


CHAMBERS AND GLOVE BOXES

BASIC GLOVE BOXES FOR MULTIPLE OPERATORS



Our catalog **#818-GBB** style Glove Box includes all the same features as the **#818-GB**, but it has greater capacity. It features a clear transparent bottom section and two (2) pairs of white Hypalon™ gloves for two (2) operators.

This is a very nice glove box for operators working in tandem (up to four [4] scientists). Refer to page #12 for product specifications.



Some optional accessory suggestions (Refer to Accessories Section).

- H.E.P.A. Filtration systems to $.3\mu$ @ 99.99%.
- Activated carbon (charcoal) filters.
- Static control/elimination systems.
- Oxygen analyzer unit.
- Germicidal (U.V.) sterilization light package.
- Additional glove ports (Up to four [4] operators).
- Computer interface ports (Blanking plates).
- Flow meters.
- Iris diaphragms for bare hand entry.
- Light sensitive units.
- Atmosphere drying set up with vacuum pump.
- Additional ground key cock valves (For inert gases).
- Cushioned transport support cart (trolley) with gas tank storage track.
- Catalyst heater unit to reduce Oxygen.
- Milled Neoprene gloves for increased sensitivity.
- Glove port plugs for damaged glove change over.
- Ten (10") inch Ø (254 mm) gloves and ports.
- The "Bright Light" illumination system (40,000 hour lamp warranty).
- Exhaust outlet for your house exhaust system.
- "VHP" sterilization requirements (With "Banjo™" fittings).
- Suspended shelf packages for increased culture dish capacity. Refer to page #38.

Did you Know?
One (1) Cubic Foot
is equal to 28.3
Liters in volume.

#818-GBB (Domestic Model)

#818-GBB/EXP (Export Model)

Refer to page #12 for specifications.

Some possible applications.

- Pharmaceutical Research Labs.
- Bio-Tech Research Facilities.
- Animal Science Research.
- Clean Room Facilities.
- Waste Management Research.
- Micro-Electronics Assembly.
- Chemical Company Labs.
- Plasma Welding Applications.
- Nano Technology Enclosures.
- Universities & Medical Schools.
- Industrial/Contract Research Applications.
- Microbiology (Anaerobic) Studies.
- Food Science Research.
- Toxicology, Biology, Geology, Zoology & Pharmacology Applications.

NOTE: Documentation for I.Q., O.Q. and F.A.T. requirements are now available from Plas-Labs. Refer to page #27. Contact us for complete details.

CHAMBERS AND GLOVE BOXES

COMPACT GLOVE BOXES FOUR SIZES AVAILABLE

The **#830-Series** of glove boxes has been developed to minimize drafts and maximize the effectiveness of your analytical balance. All units feature high visibility of the interior, cleanliness, safety, complete containment, portability, and a draft free atmosphere.

These compact glove boxes are available in four (4) user friendly sizes. They are excellent for use in hazardous situations where materials to be weighed must be contained. The chamber is ideal for determining the dry weight measure or moisture content of aqueous solutions, adhesives, cereals, toxic chemicals, pharmaceuticals, paper, plastics, and light radio active materials. All units minimize product degradation caused by moisture absorption.

They may also be used when working with toxic substances, asbestos fibers, sewage residue, and harmful liquid vapors.

NOTE: These units are not recommended for use with aromatic hydrocarbons.

When working with analytical balances, top loaders, bench scales, and/or moisture analyzers, we recommend placing those items upon a marble block for added stability.

Special access doors are available if transfer chambers are not desired. Refer to pages #15 and #17.

What's included with the glove box (Basic Features):


- Clear acrylic plastic main body (3/8" [10 mm] thick).
- Removable top section for easy installation of equipment.
- Four (4) key cock valves for purging. Two (2) on the main body and two (2) on the transfer chamber.
- All clamps adjustable to compensate for normal wear.
- White ambidextrous Hypalon™ gloves for superior chemical and U.V. resistance.
- Hospital grade multiple electrical outlet (socket) strip. UL, CSA, and CE approved.
- Small pressure relief "pop" valve minimizes glove "fight back."
- Transfer chamber is 12" (305 mm) long x 11" (280) I.D.
- Standard 24 month warranty (Not gloves or consumables).
- Small replaceable H.E.P.A. and organic vapor canister. filter. Face mask style. Refer to **#830-FILTER** on this page.

Electrical requirements (Without any accessories).

North American 115-120 Volts, 60 Hz. 12 Amps.
Export models 220-240 Volts, 50 Hz. Amps will vary as to end users country.



#830- ABC Compact Glove Box
#830- ABC/EXP (Export Model)



Small canister filter and pressure relief "Pop" valve. These are standard on all #830-Series glove boxes.

#830- FILTER (HEPA & CARBON)

#830-FILTER set up includes a small "Face Mask" style HEPA and CARBON filter canister. Sold in packages of two (2). Pop valve not included. Refer to Catalog **#830-FILTER** in Accessories Section.

CHAMBERS AND GLOVE BOXES

COMPACT GLOVE BOXES

ADDITIONAL INFORMATION

The **#830-Series** of glove boxes all include the same features and are similar in shape.

For units with access doors instead of transfer chambers, please refer to photo on [page #17](#) and [Specifications Chart](#) below.

Did you Know?
Acrylic Glove Boxes and Chambers are 12% less permeable (through the wall leakage) than Vinyl Glove boxes or Bags.



#830-ABD and #830-ABD/EXP (Export)

NOTE: All export units are shipped crated using approved & certified lumber.

Product Specifications for all #830-SERIES Units including North American and Export models.

	#830-ABC Two (2) gloves.	#830-ABB Two (2) gloves.	#830-ABD Three (3) gloves.	#830-ABE Four (4) gloves.
I.D.	28"w x 23"d x 29" h 710 mm X 580 mm x 740 mm h	35"w x 29"d x 30" h 890 mm X 740 mm x 760 mm h	48"w x 29"d x 32" h 1,220 mm X 740 mm x 810 mm h	60"w x 29"d x 32" h 1,524 mm X 740 mm x 810 mm h
O.D. Includes Transfer Chamber	43"w x 24"d x 31" h 1,100 mm X 610 mm x 790 mm h	49"w x 30"d x 31" h 1,250 mm X 760 mm x 790 mm h	63"w x 31"d x 35" h 1,600 mm X 790 mm x 860 mm h	76"w x 32"d x 35" h 1,930 mm X 812 mm x 890 mm h
Top Opening	17" x 23" 432 mm x 585 mm	24" x 29" 610 mm x 735 mm	24" x 42" 610 mm x 1,065 mm	24" x 54" 610 mm x 1,370 mm
Approximate Volume	10.8 Cubic Ft. 306 Liters	17.6 Cubic Ft. 498 Liters	25.8 Cubic Ft. 730 Liters	32.2 Cubic Ft. 912 Liters
Approximate Shipping Weight (Crated)	230 Pounds. 104 Kilos	275 Pounds. 125 Kilos	310 Pounds. 141 Kilos	400 Pounds 182 Kilos
Access Door Opening (Optional)	15.5" x 22" h 400 mm x 560 mm h	21.5" x 22" h 546 mm x 560 mm h	21.5" x 24" h 546 mm x 610 mm h	21.5" x 24" h 546 mm x 610 mm h

CHAMBERS AND GLOVE BOXES

COMPACT GLOVE BOXES WITH HEPA FILTRATION

We offer three (3) models of the **#830-SERIES** glove boxes featuring HEPA filtration packages. There are two (2) different versions available. Refer to listings 1 & 2 below. **NOTE: All three models include the same basic #830-Series features and are available in four (4) sizes. Refer also to pages #14 & 15.**

“Closed Loop Style” CLC Style.

This system continually re-circulates and filters the internal atmosphere of the glove box. It is a negative pressure containment unit. Magnehelic™ gauges indicate both inner chamber pressure and pressure drop levels across the filters. These gauges indicate when it is time to change the HEPA filters. Add 15.5” (395 mm) to the overall height of this style glove box.

What’s included with the “Closed Loop” CLC-Style glove box:

- Totally enclosed blower (fan) unit.
- Two (2) cylindrical 8” inch Ø (203 mm) HEPA filters.
- Two (2) Magnehelic gauges. Domestic models in Inches of H₂O and Export in Pascals.
- All necessary hoses and connections.
- Non-Hydroscopic pre-filters.
- Standard 24 month warranty (Not gloves or consumables).

“Open Exhaust Style” (Two filter versions)

These glove boxes draw the ambient room atmosphere through a HEPA filter and then into the chamber. Both versions are set up for negative pressure operations. Exhaust atmosphere is filtered through a HEPA or 5μ (micron) glass media filter. The “open exhaust” outlet can be set up to connect to your house exhaust system. Please specify a **#800-VENT** outlet accessory.

1. Catalog #800-HEPA/S (Single HEPA and glass media exhaust filter)
2. Catalog #800-HEPA/D (Double HEPA filters)

Variations in different types of filters are:

1. HEPA's
2. ULPA's
3. Activated Carbon (Charcoal). Four (4) different compositions. **Refer to pages #36 & 37.**



#830-ABB/CLC (Closed Loop Glove Box)
#830-ABB/CLC/EXP (Export-Closed Loop Glove Box)



1.

#830-ABC/SPCL/800-HEPA/S (Open Exhaust Style)
#830-ABC/EXP/SPCL/800-HEPA/S (Export)



2.

#830-ABC/SPCL/800-HEPA/D (Open Exhaust Style)
#830-ABC/EXP/SPCL/800-HEPA/D (Export)

CHAMBERS AND GLOVE BOXES

COMPACT GLOVE BOXES

ADDITIONAL INFORMATION AND STYLE VARIATIONS

The **#830-SERIES** of glove boxes can be modified to create different variations.

Access doors can be substituted for transfer chambers. Sizes of the doors vary according to the size of the glove box. Refer to **#1** in the photograph.

Exhaust (4" [102 mm] Ø) vents may also be added for connecting to your house exhaust system. Refer to **#2** in the photograph. Refer also to page #18 **#800-VENT/OUTLET**.

Square (box type) transfer chambers can be substituted for cylindrical (round) ones. Square transfer chambers have pressure limitations however. **NOTE:** The maximum working pressure on the square transfer chambers is 6" of WC (150 mm).

We now offer a new automatic door locking system for these transfer chambers. Catalog **#800-DOOR/ILOCK**.



#830-SHELF-ABC

Catalog **#830-SHELF-ABC**, Suspended shelf accessory. Total shelf capacity is about 150 Petri dishes depending upon chamber size. Refer also to page #38.



#830-CP/11

#830-CP/11/EXP (Export model)

Please refer to page #15 for door and glove box dimensions. To include a transfer chamber with the **#830-CP** unit order **#800-TRANS-CH**. Specify desired location of transfer chamber.



Square Transfer Chamber with #800-DOOR/ILOCK

NOTE: Maximum working pressure on this style transfer chamber is 6.0" of H₂O (0.44" of Hg).

When one door is opened, the other door locks automatically. Locking force is 120 Lbs. Refer to page #47 for more details.

CHAMBERS AND GLOVE BOXES

GERMICIDAL WORKSTATION (STILL AIR ENCLOSURE)

Our catalog **#820-TCH** Workstation has been designed to fit nicely on a lab counter top. It provides simple protection when a more sophisticated device is not required. It is an ideal enclosure for many acute and simple research procedures.

Its modest weight of 60 pounds (27.5 Kg) makes it readily transportable. The main housing is formed of one (1) piece clear acrylic and offers total viewing of your task at hand. High impact clear PETG units or amber "light sensitive" units are also available.

This Workstation can also be used as a PCR Chamber or Tissue Culture Hood.

Features included with the product:

- Standard twenty-four (24) month warranty (Not consumables).
- Formed one (1) piece clear acrylic top section .375" (950 mm) thick. With "easy clean" corners."
- "Bright Light" illumination system (40,000 hour lamp guarantee).
- Formed in place liquid barrier to facilitate easy clean up of spills.
- The U.V. system is rated at 254 nm. Automatic timer is included. Beta ray protected **but not gamma rays**.
- All controls are illuminated and the system is fused for safety.
- Ground key cock valve for introduction of Oxygen or your inert gas of choice.
- Formed plastic bottom with "spill guard." This for easy clean up.
- Hospital grade multiple electric outlet strip (socket).
- Easy access sash. Iris diaphragms included.

Electrical requirements.

North American 115-120 Volts, 60 Hz. 12 Amps.
Export Models 220-240 Volts, 50 Hz. Amps will vary as to the end users country.

#820-TCH Specifications


I.D.	41"w x 28"d x 17"h 1,040 mm x 710 mm x 432 mm
O.D. (w/Flange)	44"w x 29"d x 30"h 1,118 mm x 736 mm x 762 mm
Sash Opening	8.5" h x 27.5 long 216 mm x 700 mm
Shipping Weight (Crated)	110 Pounds 50 Kilos



#820-TCH Germicidal Workstation
#820-TCH/EXP Workstation (Export Model)



Workstation Sash for Draft Control
featuring Iris Diaphragms.



Optional Workstation Vent #800-VENT/OUTLET

Catalog **#800-VENT/OUTLET** can be mounted just about anywhere on the enclosure and is available in any diameter. Use this to connect to your house exhaust system. Standard diameter is 4"Ø (102 mm).

CHAMBERS AND GLOVE BOXES

“LA PETITE” GLOVE BOX

Plas Labs is proud to offer a new concept in small glove boxes. The **#815-PGB** “Petite Glove Box” is comprised of essentially two (2) components. The top and bottom sections are sealed using a revolutionary new gel system commonly used in the clean room industry for sealing HEPA filters.

The La Petite is only two-thirds the size of our **#818-GB** Basic Glove Box. It has been designed to be ergonomically correct for the average person.

Careful attention to research requirements and our over forty (40) years of experience helped make this glove box very user friendly and inexpensive.

Basic Features

- Lightweight and easy to move from lab to lab.
- Optically clear, one piece formed top is unbreakable.
- Bright white one piece bottom formed out of high temperature Noryl™ and has “easy clean” corners.
- Large access door (12” inch Ø [305 mm] opening).
- Vacuum gauge on transfer chamber.
- All stainless steel clamps are adjustable and have security locking capabilities.
- White leveling tray for transferring liquids.
- Modern, Gel Sealant (Refer to Gel data below).
- Large, interior Rh gauge.
- Four (4) ground key cock valves for purging and gas entry,
- Hospital grade multiple electrical outlet (socket) strip.
- Latex gloves 7”Ø (175 mm). Neoprene or Hypalon™ also available (Special order).
- Transparent and vacuum capable transfer chamber size 10”long (254 mm) x 9.75”I.D (245 mm).
- Standard twenty-four (24) month warranty (Not gloves).

Specifications

	#815-PGB
I.D.	31.5” w x 25.5” d x 21.5” h 800 mm x 650 mm x 545 mm
O.D. Includes Transfer Chamber	49” w x 29” d x 22” h 1,245 mm x 740 mm x 558 mm
Approximate volume	9.5 Cubic Ft. 269 Liters
Approximate shipping weight (Crated)	75 Pounds 35 Kilos

A Re-Engineered Product!



Catalog #815-PGB

Catalog #815-PGB/EXP (Export model)

The Gel Sealant System.

This new gel sealant was specifically developed for the clean room industry to seal HEPA filters.



Properties of the gel.

- The gel does not support biological growth.
- It will not “out gas” and is non-toxic.
- It is an air tight, water tight barrier.
- It contains no silicone that might cause contamination.
- It is highly adhesive, yet quick releasing and self healing.

NOTE: If operations are expected to exceed the pressure limitations listed below, please refer to our **#818-GB**, Basic Glove Box.

NOTE: “La Petite” Pressure Limitations.

- 2.0” WC (51 mm)
- 0.5 kPa
- 5.0 MBAR
- .077 PSI

Electrical requirements.

North American 115-120 Volts, 60 Hz. 12 Amps.
Export Models 220-240 Volts, 50 Hz. Amps will vary as to end users’ country.

CHAMBERS AND GLOVE BOXES

PCR (UV) CHAMBER

This Catalog #825-PCR(UV) has been designed to help improve the accuracy of P.C.R.[™] and general Tissue Culture procedures. The chance for air borne contamination during D.N.A. sequencing is greatly reduced.

The chamber is a “still air enclosure” which contains both our “Bright Light” illumination system and a U.V. germicidal lamp. The U.V. germicidal system is rated at 254 nm and will decontaminate all exposed surfaces in the interior. The lamp timer should be activated 15 minutes before and 15 minutes after each use.

The front viewing panel and main housing are .375” (9.5 mm) thick optically clear material. This thickness helps to prolong the service life (about 2,500 operating hours) of the unit. The front viewing panel Catalog #825-UV(FS) is attached by six (6) screws and is easily and quickly replaced.

Access to the interior of the chamber is gained by pulling swinging doors outward. New self closing hinges have eliminated the need for locks. The doors are removable, therefore do not take up valuable lab counter space. If a corner door is accidentally opened, proximity sensors turn off the U.V. system for operator safety.

Basic Features of the Product.

- “Bright Light” illumination system (40,000 hour lamp guarantee).
- Two bright white plastic (adjustable) shelves. One is tooled to store pipettors.
- Front panel is .375” thick (9.5 mm) for Beta Ray protection. **NOTE: Not Gamma Rays.**
- Removable side access doors with slip apart hinges.
- Proximity sensors on doors for operator safety.
- Automatic timer to activate U.V. sterilization procedures.
- Bottom tray has a formed in place “spill guard” for easy cleaning.
- The U.V. system is rated at 254 nm.
- Standard 24 month warranty (Not gloves or consumables).
- Side and back walls are one piece formed optically clear acrylic .375” thick (9.5 mm).
- Main housing and top are removable for installation of large pieces of equipment.



#825-PCR(UV) Chamber &
#825-PCR(UV)EXP Chamber

Product Specifications

I.D.	23.5”w x 17”d x 21” h 597 mm x 432 mm x 530 mm
O.D.	24”w x 18”d x 28”h 610 mm x 457 mm x 711 mm
Shipping Weight (Crated)	110 Pounds 50 Kilos

NOTE: If HEPA filtration is preferred, refer to Catalog #825-PCR/HEPA illustrated on page #21.

Electrical requirements.

North American 115-120 Volts, 60 Hz. 12 Amps.
Export Models 220-240 Volts, 50 Hz. 6 Amps.

CHAMBERS AND GLOVE BOXES

HEPA FILTERED PCR CHAMBER

NEW PRODUCT !!



The new HEPA filtered PCR Chamber, Catalog #825-PCR/HEPA, includes all of the same basic features as our standard #825-PCR(UV) Chamber, but includes a factory installed HEPA filter system. The HEPA filter system is rated at 99.9997% efficiency at 0.3 (μ) microns. No tools are required for HEPA filter change over.

This is a positive pressure research chamber. The fan (blower) noise level is rated at less than 40 dBA.



#825-PCR/HEPA
NOTE: HEPA housing on rear wall.

#825-PCR/HEPA
#825-PCR/HEPA/EXP (Export Model)

Product Specifications

I.D.	23.5" w x 17" d x 21" h 597 mm x 432 mm x 530 mm
O.D.	26" w x 23" d x 28" h 610 mm x 457 mm x 711 mm
Shipping Weight (Crated)	112 Pounds 51 Kilos

Additional Features of the Basic #825-PCR. (Refer also to page #20.)

- "Bright Light" illumination system (40,000 hour lamp guarantee).
- "Mini HEPA" filter size 4" x 9.5" x 1" thick (101 mm x 240 mm x 25.4 mm thick) Replacement part #MS2199.
- No tools required to change over the HEPA filter.
- Fan (ventilator) noise level is rated at less than 40 dBA.
- Formed in place liquid barrier to facilitate easy clean up of spills.
- Standard twenty-four (24) month warranty (Not consumables-filter).

Electrical requirements.

North American 115-120 Volts, 60 Hz. 12 Amps.
Export Models 220-240 Volts, 50 Hz. 6 Amps.

CHAMBERS AND GLOVE BOXES

LOW VELOCITY FUME HOODS (IMPROVED FEATURES)

ASHRAE 110 Compliant

This Low Velocity Fume Hood system has been developed in response to customer requests and suggestions for improvements to the original product.

These units can be used with your in house exhaust system or our accessory **#900-EXTRACT**, Fume Extraction Unit. All units are ASHRAE 110 Certified. You can select either HEPA or Activated Carbon filters (or a combination of both) for the Extraction Unit. Three (3) styles of Carbon are available.

Refer to page #40.

Our hoods are unique in that they include both horizontal and vertical “air foils” taking advantage of a “venturi” effect. They also feature rounded rear corners to guide and facilitate inner atmosphere exhaust efficiency.

The exhaust hose has a “smooth bore” interior, which minimizes air flow exhaust restrictions, is chemically resistant, FDA approved, and is included with the main work housing. It may be attached directly to your house exhaust system.

All units include the following:

- Standard 24 month warranty (Not consumables-filters).
- Units are fabricated out of .375” (9.5 mm) thick optically clear acrylic plastic.
- Hinged (two [2] position) front sash with draft deflector.
- Formed horizontal lower front airfoil.
- Two (2) formed vertical air foils with draft deflectors.
- Bright white acrylic bottom work surface.
- Two (2) side wall cable ports at 1.5”Ø (38 mm).
- Six (6”) inch Ø (about 152 mm) exhaust hook up.
- The exhaust hose is smooth bore inside and is six (6”) inch (152 mm) Ø x eight (8’) feet (2,440 mm) long.

Did you Know?

The air we breathe is only 20.9% Oxygen. The rest of our air is made up of various other gases such as Nitrogen.



#900-LVFH with #900-EXTRACT and #900-CART-LVFH

Dimensions of the LVFH units:

#900-LVFH/24

Inner work surface @ 23.250”w (590 mm) x 19”d (483 mm)
Outer dimensions @ 24”w (610 mm) x 27”d (685 mm) x 28”h (712 mm)

#900-LVFH/36

Inner work surface @ 35.250”w (895 mm) x 19”d (483 mm)
Outer dimensions @ 36”w (915 mm) x 27”d (685 mm) x 28”h (712 mm)

#900-LVFH/48

Inner work surface @ 47.250”w (1,200 mm) x 19”d (483 mm)
Outer dimensions @ 48”w (1,219 mm) x 27”d (685 mm) x 28”h (712 mm)

#900-CART-LVFH

All carts (trolleys) are the same height @ 32” (813 mm) and the same depth @ 27” (685 mm).

NOTE: The cart (trolley) width will vary according to the same width as the Fume Hood.

Electrical requirements for all components:


North American 115-120 Volts, 60 Hz. 15 Amps.
Export models 220-240 Volts, 50 Hz. 7.5 Amps.

CHAMBERS AND GLOVE BOXES

LOW VELOCITY FUME HOOD COMPONENTS

Optional types of Fume Extraction filters. Specify preferred type when ordering:

- **#900-LVFH/HEPA/CARB**, Combination HEPA and Non-Impregnated Carbon. Absorbs most aliphatic and aromatic hydrocarbons, solvents, organic vapors, aldehydes, ketones, alcohols, organic acids, esters, sulfur compounds, and odors.
- **#900-LVFH/CARB (A)**. Impregnated Carbon filter for alkaline type fumes.
- **#900-LVFH/CARB (B)**. Impregnated Carbon filter for ammonia or amines.
- **#900-LVFH/CARB (C)**. Non-Impregnated Carbon filter for absorption of aromatic hydrocarbons, organic vapors, ketones, alcohols, organic acids, esters, Nitrogen compounds, and odors.
- **#900-LVFH/HEPA**. The main HEPA (particle) filter for the Extraction Unit. Rated at 99.99% at 0.3 μ (microns).



#900-LVFH
Close up of Air Foils

NOTE: Air foils are:

- #1. Lower front horizontal deflector.
- #2. Vertical air foils (2).
- #3. Top sash draft deflector.

All ASHRAE 110 Compliant



#900-EXTRACT
Portable fume extraction unit (Optional)

The Fume Extraction Unit is a variable speed suction housing that can include your choice of filters such as HEPA or Impregnated Carbon. Filter sizes are 12" (305 mm) x 24" (610 mm).

Suction speeds can be adjusted from "0" up to 250 CFM face velocity. As the filter resistance builds up, you can compensate for it by increasing the suction speed.



#900-CART-LIFT
Adjustable Electric Lift Tables (Optional)

The #900-CART-LIFT Electronic Height Adjustment Table has been designed to provide ergonomically correct operating heights to ensure maximum operator comfort. Widths will vary according to the width of the fume hood selected. Refer to page # 22.

CHAMBERS AND GLOVE BOXES

ISOLATION AND CONTAINMENT UNITS

Isolation Systems Engineered for Specific Research Projects.

Select and start with a basic housing size and add your desired features to develop the exact isolator for your research needs.

Our rigid acrylic isolators provide complete transparency (no hidden corners), long term strength (durability), and freedom from through the wall leakage of atmosphere. Permeation rate is minimized. Resistance to permeation is 12% better than flexible vinyl style units.

The optical quality and durability of acrylic does not become brittle and crack with age. PETG units are also available as an option if VHP sterilization is anticipated. “Banjo™” type fittings are included for VHP applications.

Multiple station units offer the benefit for greater productivity and efficiency during sequential procedures. Units may also be attached in train if required. Refer to Catalog #800-T Transfer Chamber “T” in the Accessories Section. [Page #48.](#)

Features included with each isolator:

- Standard 24 month warranty (Not gloves or consumables-filters).
- The top and bottom sections are each formed of one (1) piece optically clear acrylic .375” (950 mm) thick.
- Joining gaskets are closed cell neoprene and doubled layered 1” (25.4 mm) thick for long term durability.
- Stainless Steel quick acting (adjustable) clamps and strikes.
- Ambidextrous, white Hypalon™ gloves for superior U.V. and chemical resistance.
- Stainless Steel glove port rings for added security.
- Transfer chamber at 12” long (305 mm) x 11” Ø (280 mm) I.D. NOTE: Other sizes available.
- Two (2) ground key cock valves for purging the transfer chamber.
- White leveling tray for transferring liquids into the main chamber.
- Rigid doors on the transfer chamber allow for complete purging and include adjustable stainless steel clamps.
- Rounded “Easy Clean” corners on the top and bottom sections.

NOTE: Dimensions and product details for all Series #800-RPL isolators are listed on [page #25.](#)



#805-RPL Isolator (Double Rack System)

Electrical requirements:

Amp requirements will vary as to the features selected.

North American 115-120 Volts, 60 Hz.

Export Models 220-240 Volts, 50 Hz.

The entire isolator system can be safely relocated from lab to lab by merely unplugging the electrical connection.



Catalog #800-SHELF-II (Sliding shelf accessory)

The sliding shelf accessory greatly increases the holding capacity of the isolator. Our shelf is designed to keep the valuable floor area uncluttered. No shelf supports eliminate cleaning problems on the floor. Shelf can be pushed to rear of the isolator when not in use. Refer to Accessories Section for additional shelf details. [Page #38.](#)

CHAMBERS AND GLOVE BOXES

ISOLATION UNITS & STERILITY ISOLATOR DATA

Our Sterility Isolator, Catalog #881-SI has been designed to provide a sterile, particle free working environment. The HEPA filtration system maintains a slight positive pressure in the main chamber ensuring complete product isolation. The HEPA fan unit has an adjustable speed control.

The isolator can be sterilized using most of the commercially available sterilants. If needed, VHP set ups are also available. **NOTE:** Port plugs must be used for VHP procedures.

Access to the main working chamber is achieved through our Hypalon™ glove system. Our unique bare hand entry port system is also available if preferred.

All isolator systems can be safely moved from lab to lab by merely unplugging the electrical connection.

Sterility Isolator (#881-SI) Features Include:


- Standard 24 month warranty (Not entry system or consumables).
- “Bright White” bottom work surface .375” thick (950 mm).
- Large HEPA filter size 24” (610 mm) x 48” (1,220 mm) x 6” (152 mm) thick. Same size as inner work surface area.
- Optically clear .375” thick (950 mm) acrylic construction.
- Large 10”Ø (254 mm) bare hand entry system.
- Large square transfer chamber with two (2) purging valves.
- Multiple electric outlet (socket) strip, hospital grade.
- Removable shelf.
- Adjustable speed control (70-140 FPM).
- “Bright Light” system (40,000 hour lamp guarantee).
- U.V. sterilization system (254 nm).

Some available accessories:

- Iris diaphragms.
- Flow meter set ups (PSI & Pascals) [Refer to page #46.](#)
- Static elimination Ionizers (Draft free).
- “Banjo™” type fittings for VHP sterilization.
- Glove port plugs. [Refer to page #39.](#)
- “Dip (Dunk) Tank” entry system with lid.
- Automatic Transfer Chamber “interlock” system. [Refer to page #47.](#)
- Stainless steel cushioned cart (Trolley).

Available Certification services. [Refer to page #27.](#)

1. I.Q. & O.Q. Certification.
2. Factory Acceptance Test (F.A.T.) Certification.
3. HEPA Filtration Certification (By outside agency).



#881-SI, Sterility Isolator

#881-SI/EXP, Export Model Sterility Isolator

Electrical requirements:

North American 115-120 Volts, 60 Hz. 10 Amps.
Export Models 220-240 Volts, 50 Hz. 5 Amps.

Product Specifications

Model Number	Approximate Volume	Isolator Dimensions I.D. Base and Height	Approximate Shipping Wt.
#881-SI	30.2 Cu. Ft 854 Liters	58" w x 29" d x 31" h 1,472 x 737 x 787 mm Overall Dimensions (O.D.) 78" w x 30" d x 45" h 1,982 x 762 x 1,144 mm	180 Pounds 83 Kilos
#800-RPL	20.7 Cu. Ft 586 Liters	30" w x 30" d x 35" h 760 x 760 x 790 mm Overall Dimensions (O.D.) 40" w x 40" d x 36" h 1,010 x 1,010 x 910 mm	145 Pounds 67 Kilos
#802-RPL	36.5 Cu. Ft 1,033 Liters	48" w x 32" d x 34" h 1,210 x 810 x 840 mm Overall Dimensions (O.D.) 58" w x 42" d x 36" h 1,470 x 1,070 x 915mm	186 Pounds 85 Kilos
#805-RPL	40 Cu. Ft 1,132 Liters	56" w x 32" d x 35" h 1,420 x 810 x 840 mm Overall Dimensions (O.D.) 64" w x 43" d x 36" h 1,620 x 1,090 x 915 mm	197 Pounds 90 Kilos
#806-RPL	48 Cu. Ft 1,358 Liters	72" w x 32" d x 35" h 1,830 x 810 x 860 mm Overall Dimensions (O.D.) 81" w x 41" d x 36" h 2,060 x 1,040 x 915 mm	235 Pounds 106 Kilos

CHAMBERS AND GLOVE BOXES

BIOTOXIN CONTAINMENT & DECONTAMINATION CHAMBERS

Some Samples of Custom Built Isolators.


Illustrated are two special units coupled together for biotoxin containment (left) and a decontamination chamber. They are negative pressure units and include three different special filter set ups.

- #1. Molecular Sieve.
- #2. Activated Carbon.
- #3. HEPA Filters.

Isolators are available with interlocking stainless steel support carts (trolleys) with Marine Edge tops. These tops are a safety feature in case of liquid spills.

Units as shown are set up for connection to house exhaust systems.

Refer to pages 12 & 15 for dimensions.



#818-GBB/SP (left) & #830-ABB/SP (right)

What can we build for you?

DRUG PROCESSING STATIONS & ISOLATORS SET UP IN TRAIN

Illustrated are two special #818-GBB glove boxes coupled together used as a drug processing station for multiple operators.

Units are set up for connection to house exhaust systems.

Isolators are available with interlocking stainless steel support carts (trolleys) with Marine Edge tops or safety.

Refer to pages 12 & 15 for dimensions.



#818-GBB/SP (double set up)

Note: All isolators can be constructed of clear PVC or PETG if VHP sterilization is planned. Banjo™ fittings are included. Refer to page #46.

Available Certification services.

- 1. I.Q. & O.Q. Certification.
- 2. Factory Acceptance Test (F.A.T.) Certification.
- 3. HEPA Filtration Certification (By outside agency).

CHAMBERS AND GLOVE BOXES

CUSTOM CAPABILITIES AND PROJECTS

Some Additional Samples of Custom Built Isolators.

PLAS ■ LABS is committed to providing high quality research equipment at a fair price. We realize that unique research protocols may require individualized equipment. We offer the following special services:

- Design Engineering.
- New Product Development.
- I.Q. & O.Q. Documentation.
- Custom Service/Maintenance Schedules.
- Customized Operator Manuals.
- Third Party Certification.
- CAD Drawings.
- F.A.T. Documentation.
- On Site Training.
- Various Materials of Construction, Eg., Light Sensitive Protection, Stainless Steel, Static-Dissipative Acrylic, Clear PVC, and PETG.

Our design engineering team adheres to the most up to date guidelines from the following organizations:

- International Society for Pharmaceutical Engineering.
- American Glove Box Society.
- American Biological Safety Association.
- U.S. National Institutes of Health.
- American Association for Lab Animal Science.
- U.S. Centers for Disease Control.
- The World Health Organization (The Laboratory Biosafety Manual).



Conditioning and Weighing Diesel Particulate Filters

Two stainless steel glove boxes with automatic temperature and dew point controls, blower, and HEPA filter system for auto emission studies.



A Special "Nanotechnology" Unit

Compact glove box with temperature and humidity control with data logging capabilities. Designed for use in Nanotechnology studies.



Small Tablet Press Isolator

Isolator designed to house a small tablet press. Includes inlet/outlet HEPA system, hermetically sealed feed throughs, removable top and transfer chamber



Penicillin Processing Isolator



A Drug Transfer Isolator

Isolator for transferring pharmaceutical tablets from a press to a table coater. Includes sanitary fittings, spray ball with quick disconnect and transfer chamber.

CHAMBERS AND GLOVE BOXES

ACRYLIC DESICCATOR CABINETS

Clear and Light Sensitive Units.

Our Desiccators are designed to protect, store, or assist in transporting critical materials and prevent their contamination.

They are dust and moisture free cabinets that safeguard valuable research specimens or sensitive electronic components. You may easily observe the status of your components without opening the door.

Under a slight positive N₂ charge they prevent moisture from penetrating the molecular structure of stored components.

Amber colored units “CGA” are very effective for protection from harmful U.V. rays for light sensitive components.

NOTE: #861-CG and #861-CGA doors open from the side. All other models open from the top and bottom.

All Desiccators Include These Features

- Units are all stackable with non-slip “grip tight” feet.
- Standardized construction means interchangeable shelves.
- All shelves have perforated slots (holes) for optimum gas saturation.
- All Desiccators with two or more cubicles have a built in divider for added strength.
- Units are constructed of optically clear “jeweled” acrylic.
- Slip-a-part hinges are standard.
- Sealing clamps are springloaded for uniform locking tension.
- Built in hygrometer (Rh) gauge.
- All units have “door guards” to prevent warping and possible loss of internal atmosphere.
- Closed cell “Poron” gasketing does not take a set.
- All shelves spaced at 2” (50.8 mm).
- Special “gasket guard system” prevents gaskets from taking a set and causing a leak.
- Each cubicle includes two (2) Nylon gas serrated nipples 1/8” & 1/4” (90 degree elbows with hose barbs).



#860-CGA, #860.5-CG, #861-CG, #862-CG

Static dissipative or clear PVC units are also available.



#800-FLOW/D

Useful for careful adjustment of flow rate.



#800-PRV (Pressure Relief)

Minimizes the chance for explosion or implosion.

Refer to [page #46](#) for these accessories specifications.

Desiccator Specifications

Model Number	Color	Number of Cubicles	Number of Shelves	Number of Shelf Positions	Desiccator Sizes W x L x H
860.5-CG 860.5-CGA	Clear Amber	One (1)	Two (2)	Two (2)	12 x 12 x 6" 305 x 305 x 152 mm
860-CG 860-CGA	Clear Amber	One (1)	Two (2)	Five (5)	12 x 12 x 12" 305 x 305 x 305 mm
861-CG 861-CGA	Clear Amber	Two (2)	Four (4)	Ten (10)	12 x 12 x 24" 305 x 305 x 610 mm
862-CG 862-CGA	Clear Amber	Two (2)	Four (4)	Ten (10)	18 x 18 x 18" 457 x 457 x 457 mm
863-CG 863-CGA	Clear Amber	Two (2)	Four (4)	Ten (10)	24 x 18 x 24" 610 x 457 x 610 mm

CHAMBERS AND GLOVE BOXES

MULTIPLE CUBICLE DESICCATOR CABINETS

These multiple cubicle desiccator storage cabinets are much safer than individually stacked units. All of these storage units are mounted on stainless steel support carts (trolleys) with casters and remain stable during transportation.

Transparent amber colored units provide U.V. protection for light sensitive materials.

Three (3) standard sizes are available and custom fabrication is welcomed.

Benefits of the storage cabinets.

- Complete visibility of stored contents.
- Unitized construction eliminates “tip overs”.
- Safe air tight storage of critical components.
- U.V. degradation is eliminated using amber units.
- Easy to transport from lab to lab.
- Ideal for Clean Room applications.
- Moisture sensitivity problems are eliminated.

What's included with desiccator storage cabinets.

- Standard 24 month warranty.
- Clear or amber models.
- Unitized construction for greater durability and strength.
- Three standard cubicle set ups (sizes) are available.
- Two (2) shelves with slots included per cubicle
- Two (2) Nylon gas nipples included per cubicle.
- Stainless steel support frames (trolleys) with casters, brakes and locks.
- Interior volume of each individual cubicle is 3.8 Cu. Ft (107 L).
- Closed cellular (one piece) Poron gasket.
- All doors have “anti-warp” framework and heavy duty “swell” latches.
- Exterior rear plenum with flow meter.

Optional accessories (Refer to Accessories Section).

- Controlled Purge System, Catalog #800-CPS.
- Pressure relief valve, Catalog #800-PRV.
- Additional shelves.
- “Anti-Sag” supports allow for extra weight capacity.
- Extra desiccant storage trays.
- Small or large flow meters.

NOTE: Interior dimensions of all cubicles are 23.5” x 23.5” x 11.75” h (596 mm x 596 mm x 298 mm h).

All door openings are 22.50” x 9.25” h (570 mm x 235 mm h).



#860-CGA/10 (For light sensitive materials)



#860-CG/8 (Clear unit)

CHAMBERS AND GLOVE BOXES

LAB (CO₂) CHAMBERS, HEAVY DUTY UNITS

The modest cost of these units enables the investigator to easily dedicate one chamber for one specific research project.

These superb cabinets have superior strength and sealing characteristics, utilizing heavier gauge transparent acrylic plastic with formed top corners to minimize section joints.

The structural integrity of the cabinets ensures economy in the cost of operation. After the initial gas purge to clear a cabinet loaded with sensitive components, gas consumption is reduced to an absolute minimum. This is because the cabinet will not flex or deform resulting in gas leakage. When pressurized with an inert gas, no significant change in pressure level occurred over a five (5) day test.

Lab Chamber Features:

- Standard 24 month warranty.
- Heavy gauge .375" and .500" thick transparent jeweled acrylic plastic.
- Spring steel tension clamps for uniform gasket pressure.
- Non slip feet for stacking.
- Closed cellular (one piece) Poron gasket.
- Unique "gasket guard" system to prevent gaskets from taking a "set."
- Lightweight, easily transportable. Handles are available.
- Shelves are spaced at 2" (50.8 mm).
- Adjustable perforated shelves for complete atmosphere movement throughout the chamber.
- Nylon inlet and outlet gas nipples, hose barb style, size 1/8" & 1/4" hose barbs.



#800-FLOW/D (Accessory)

Optional Flow Meters enable the operator to establish and monitor a constant flow of CO₂ at lower volumes.

Refer also to page #46.



#850-LCS and #850-LCM

Lab Chamber Specifications

Catalog Number	Sizes W x L x H	Number of Shelves	Number of Cubicles
850-LCS	12" x 12" x 6" 305 mm x 305 mm x 152 mm	Two (2)	One (1)
850-LCM	12" x 12" x 12" 305 mm x 305 mm x 305 mm	Two (2)	One (1)
850-LCL	12" x 12" x 24" 305 mm x 305 mm x 610 mm	Four (4)	Two (2)
850-1L	18" x 18" x 18" 457 mm x 457 mm x 457 mm	Two (2)	One (1)
850-2L	24" x 18" x 24" 610 mm x 457 mm x 610 mm	Four (4)	Two (2)
850-3L	24" x 24" x 24" 610 mm x 610 mm x 610 mm	Four (4)	Two (2)
850-4L	36" x 36" x 36" 914 mm x 914 mm x 914 mm	Eight (8)	Four (4)

NOTE: Amber plastic is available to protect your light sensitive materials.

Clear PVC or static dissipative plastic also available.

CHAMBERS AND GLOVE BOXES

STAINLESS STEEL DESICCATORS

We are pleased to offer various sizes of heavy duty stainless steel desiccator cabinets. These units are constructed of highly polished #304 stainless steel and feature unitized construction.

We have over forty (40) years experience building thermoplastic desiccators. This knowledge is now put to use in our stainless steel units.

They are ideal when gathering samples in the field, or just moving materials safely from lab to lab. They are unbreakable.

We welcome custom sizes and you may choose either safety glass or static dissipative plastic viewing panels on the doors.

Sensitive Flow Meters and other accessories such as grounding features are available.

Stainless Steel units include:

- Standard 24 month warranty.
- Constructed of 16 Gauge, 304 stainless steel.
- Full length stainless steel hinge.
- Adjustable stainless steel sealing clamps.
- Safety glass front viewing panel.
- Adjustable shelf positions.
- Two (2) stainless steel shelves included.
- Two (2) ground key cock gas valves.
- Non skid foot pads.
- Closed cellular silicone gasketing.
- Optional formed white plastic desiccant tray.
- Flow meter (0-10 LPM) included.
- Pressure relief valve (Catalog #800-PRV) included.

Benefits of the Stainless Steel units:

- Units are easily grounded if necessary.
- Stainless Steel construction eliminates breakage.
- Increased service life (durability).
- Higher levels of vacuum/pressure are attainable.
- Increased chemical resistance.
- "Outgassing" is eliminated.

Did you Know?
In the U.S.A. there
are about 2,080 work
hours in one
(1) year. (Less
coffee breaks!)



#861-SS with optional #800-FLOW/D with gas valves.

Refer to page #46 for Flow Meter Specifications.

Stainless Steel Desiccator Specifications

Model Number	Number of Shelves	Interior Size W x D x H	Approximate Shipping Wt.
860-SS	2	12" x 12" x 12"h 305 mm x 305 mm x 305 mm h	40 Pounds 18.5 Kilos
861-SS	2	12" x 18" x 24"h 305 mm x 457 mm x 610 mm h	50 Pounds 23 Kilos
862-SS	2	20" x 18" x 30"h 508 mm x 457 mm x 762 mm h	85 Pounds 40 Kilos



#800-FLOW/D

Refer to page #46 for specifications.

CHAMBERS AND GLOVE BOXES

“INTEGRITY” STAINLESS STEEL GLOVE BOX

PLAS ■ LABS is pleased to offer our re-engineered stainless steel glove box which is based upon our over forty (40) years of experience. This rugged unit is ideal for applications in the field as well as the lab. Catalog #819-SS/GB.

Unitized construction provides complete integrity for protection of operators. Samples can be gathered and studied in the field or moved back to your laboratory as needed. Complete integrity is maintained.

Product benefits:

- Complete portability from lab to lab.
- Easy to load large bulky equipment.
- The unit does not “out gas”.
- Increased chemical resistance of stainless steel.
- Excellent for low humidity (Rh) requirements.
- Neutral (non-reactive) environmental work.
- Complete safety for transferring anaerobic field samples.
- Units can be coupled together for multiple user operations.

Two sizes of Transfer Chambers are available.

1. Standard 11” I.D. (280 mm) x 12” long (305 mm)
2. Special size 17” I.D. (455 mm) x 12” long (305 mm)

Electrical requirements:

North American 115-120 Volts, 60 Hz. 15 Amps.
Export Models 220-240 Volts, 50 Hz. 7.5 Amps.

A Re-Engineered Product !



Catalog #819-SS/GB

Catalog #819-SS/GB/EXP (Export model)

Product Data	Catalog #819-SS/GB
I. D.	40”w x 29”d x 27” h 1,020 mm x 740 mm x 690 mm h
O.D. Includes Transfer Chamber	53”w x 30”d x 34” h 1,350 mm x 765 mm x 864 mm h
Approximate Volume	17.3 Cubic Ft. 489 Liters
Approximate Shipping Weight (Crated)	275 Pounds 125 Kilos



Catalog #819-SS/GB Interior View

What’s included in the glove box:

- Standard 24 month Warranty (Not gloves or consumables).
- Rugged 12 Ga. stainless steel housing.
- Safety glass viewing port.
- Optically clear transfer chamber for viewing samples.
- Four (4) ground key cock valves for purging.
- Total front panel access for large equipment.
- Adjustable stainless steel hardware throughout.
- Ambidextrous white Hypalon™ gloves (10”Ø [254mm])
- “Bright Light” system (40,000 hour lamp guarantee).
- Polycarbonate inner transfer chamber door. (No-out-gassing).
- Multiple electric outlet (socket) strip, hospital grade.
- Dual (gas piston) lifters on the front door.
- Neoprene gasketing standard.

CHAMBERS AND GLOVE BOXES

MULTI-PURPOSE VACUUM CHAMBERS

These fine vacuum chambers are similar in design to our glove box transfer chambers that have been used for over forty (40) years.

Custom sizes and shapes are available.

Product benefits include:

- Low to medium vacuum maintenance. Down to -25" of Hg.
- Useful for degassing and freeze drying processes.
- Easy for transporting contamination and moisture sensitive components.
- Can be used in the vertical or horizontal positions.
- Complete transparency gives visual awareness of your project materials.
- Eliminates the need messy vacuum grease.

Vacuum Chamber features include:

- Standard 24 month warranty.
- Two (2) ground key cock valves for purging.
- Main chamber is removable from the white base.
- Formed "bell shaped" ends for greater safety.
- Heavy duty .500" thick jeweled acrylic.
- White base includes non-slip feet.
- Vacuum gauge.
- Greaseless Neoprene "O" ring.
- Stainless steel hardware.
- Bright white leveling tray.

Useful Accessories:

- **#800-PUMP**, Diaphragm Vacuum-Pressure pump (110-or 220/240 volts).
- **#800-PRV**, Pressure Relief Valve. Lessens chance for implosion or explosion.



#800-VAC
Horizontal Configuration-Set Up



#800-VAC
Vertical Configuration-Set Up

Catalog Number	Approximate Volume	Size
#800-VAC/S	1,350 Cubic In. 22 Cu. Liters	14" x 12"Ø 355mm x 305mm
#800-VAC/L	3,200 Cubic In. 52.5 Cu. Liters	14" x 18"Ø 355mm x 457mm

Accessories for all Chambers and Glove Boxes

MAJOR DEVICES

PLAS ■ LABS offers a complete line of beneficial accessories for chambers, isolators, and glove boxes. Select the important accessory that will help you in your research success!

"We produce research equipment that helps you produce results!"



#800-DOI, Digital Oxygen Indicator
continuously displays and monitors chamber PPM levels.

Accuracy: $\pm 1\%$ of full scale.
Resolution: 0.1ppm for ranges <100 ppm.
1.0 ppm for ranges >100 ppm.

Electrical requirements for all products this page:

North American: 115-120 Volts
60 Hz. 5 Amps.

International (Export): 220-240
Volts 50 Hz. 2.5 Amps.

#800-DOI, Digital Oxygen Indicator
#800-DOI/EXP Export Model




#800-AS/SPI "Counter Spi" Ionizer
#800-AS/SPI, Export Model

Emits alternating positive and negative Ions. This is a non-air assisted static eliminator (no blower-fan). No particle contaminant's are introduced into the chamber or glove box. This is a perfect device for working with fine powders and chemicals. Effective range is 36" (915 mm).

#800-HCB, Humidity Control Box
#800-HCB, Export Model

This device accurately controls and maintains Rh levels inside any our of chambers, isolators, or glove boxes.

NOTE: Order #800-HCB/IG for inert gas dehumidification applications.



#800-HCB/IG, Humidity Control Box
Rear View

Accessories for all Chambers and Glove Boxes

MAJOR DEVICES (CONTINUED)



Major NEW products !!

#800-CPS, Controlled Purge System
#800-CPS/EXP Export Model


This device has been designed to monitor and control positive pressure within any of our glove boxes, isolators, or desiccators.

The unit features a programmable logic controller that continuously displays the chamber pressure as well as the desired "set point." Once you have entered your desired set point, (from 0.0 to +2.0" WC or 0 to +500 Pa), the controller will automatically maintain that point.



#800-CPS/O2 and #800-CPS/O2/EXP

This accessory is specifically for Oxygen use only.
It includes all the features of the basic unit above.



#800-CPS/O2 and #800-CPS/O2/EXP (rear view)




#800-BACTI (Microcinerator™)
#800-BACTI/EXP Export Model

The BACTI Microcinerator™ sterilizes metal loops and needles without the need for an open flame.
Aerosaling of micro-organisms is eliminated.

The Microcinerator™ protects technicians from dangerous gases, flames and spatter. Loops or needles are sterilized within 5-7 seconds.

The ergonomic design reduces hand fatigue and allows the technician to access the burner from five (5) various positions. Suction cup feet provide great stability!



#800-BLS "Bright Light" System
#800-BLS/EXP Export Model

Our new "Bright Light" system is guaranteed to have a 40,000 hour bulb life. Illumination of the glove box or isolator is dramatically increased. Heat production is minimized making the unit much more comfortable in which to work.

Electrical requirements for all products this page:

North American: 115-120 Volts 60 Hz. 5 Amps.

International (Export): 220-240 Volts 50 Hz. 2.5 Amps.

Accessories for all Chambers and Glove Boxes

HEPA FILTRATION VARIATIONS & COMPONENTS

We offer HEPA filtration setups in several different configurations depending upon your application need. All setups include our **#800-MAGNA**, Magnehelic Gauge indicators.

HEPA filters can be in the following shapes:

- A. **Round. Catalog #800-H8.** Our standard “Multi-Purpose” filter.
- B. **Rectangular.** For PCR and special isolators.
- C. **Canister (Capsule).** Catalog #800-HEPA/P.
For filtering pressure overflows.

NOTE: Our Carbon (Charcoal) filters can also be mounted onto our HEPA filters. Four (4) styles available. Refer to photo this page.



#800-HEPA/P Canister (Capsule style)


Can be used for chamber pressure overflow and/or mounted on our diaphragm pump.

Traps impurities before they can be exhausted into the room atmosphere.




#800-HEPA/P mounted on #800-PUMP
For export use specify #800-PUMP/EXP

Electrical requirements for all products this page:
North American: 115-120 Volts 60 Hz. 5 Amps.
International (Export): 220-240 Volts 50 Hz. 2.5 Amps.




#800-MAGNA, Magnehelic Gauge-In Pascals (Export)
#800-MAGNA, Magnehelic Gauge-In Inches of H₂O

These gauges are used to monitor filter effectiveness and efficiency. They indicate when it is time to change your filter(s). Refer also to page #46 for more details on operation.



Typical Carbon(Charcoal) filter. Four (4) Carbon compositions available. Refer to pages #37 & 40.

#800-H8 Typical round HEPA (replacement) filter. Eight (8") Ø (203 mm)



#800-HEPA/S (Single filter typical set up)

The single filter set up includes 8"Ø (203 mm) HEPA filter, .4 µ micron pre-filter, blower (fan), mounting plates, and **#800-MAGNA set** (Refer also to page #37 & #46.)

Accessories for all Chambers and Glove Boxes

HEPA FILTRATION VARIATIONS & COMPONENTS (CONTINUED)



#800-HEPA/S and #800-HEPA/S/EXP for export models (External view with exhaust outlet connector)



#800-HEPA/S with filter access door open
Internal view (When changing filters)

Additional Filtration components (not shown).

#CF-50/9 Carbon Filter “Spun Media” pad (9”Ø) Color black. Twelve (12) per package.

#GF-50-B9 Glass media Filter Pad (9”Ø) Color pink, five (5 μ) microns. Twelve (12) per package.

#MF-50/9 Spun Polypropylene Pre-filter Pad (8”Ø) diameter. Color white, point four (.4 μ) microns. Twelve (12) per package.



#800-HEPA/D. Open exhaust (double)
HEPA filtration system


NOTE: Refer also to page #16 for more details on double HEPA filtration systems.

Special Carbon filter types available.

#800-CARBON/A Replacement Carbon Filter (8”Ø) (203 mm) diameter. One (1) per package. Impregnated Carbon Filter for alkaline type fumes.

#800-CARBON/B Replacement Carbon Filter (8”Ø) (203 mm) diameter. One (1) per package. Impregnated Carbon Filter for ammonia or amines.

#800-CARBON/C Replacement Carbon Filter (8”Ø) (203 mm) diameter. One (1) per package. Non-Impregnated Carbon Filter for absorption of aromatic hydrocarbons, organic vapors, keytones, alcohols, organic acids, esters, Nitrogen compounds, and odors.



#800-HEPA/S Typical single filter set
up on the #830-SERIES glove boxes

Accessories for all Chambers and Glove Boxes


SUSPENDED SHELF PACKAGES

Suspended (sliding) Shelf packages can be added to most of our various glove boxes, isolators and chambers. They greatly increase the culture plate holding capacity of the chamber.

Sliding shelves can be pushed to the rear of the chamber when not in use.

#800-SHELF-I, Suspended Sliding Shelf for Single Operator units. Plate capacity is about 340. Fabricated of .375 thick clear acrylic, with shelves spaced at 7" (178 mm) apart. The entire assembly is suspended from the chamber top to keep the floor area uncluttered. Overall size is 18" (457 mm) high x 18" (457 mm) wide x 10" (254 mm) deep.

#800-SHELF-II, Suspended Sliding Shelf for Multiple Operator units. Plate capacity depends upon model used. Fabricated of .375 thick clear acrylic, with shelf spacing 9" (228 mm) apart. The entire assembly is suspended from the chamber top to keep the floor area uncluttered. Overall size is 22" high (560 mm) x 16" deep (406 mm). Length depends upon chamber model selected.



#800-SHELF-II (Multiple Station Glove Box)

NEW!



#800-SHELF-830 (Shelf for #830-SERIES units)

#800-SHELF-I, For chamber models #818-GB, #850-NB, #855-AC, #888-SIMPLICITY, and #890-THC.

Plate Holding Capacity if dish size is 100 x 15 mm.

Top tier capacity (5 plates high).	=	60 dishes
Tiers 2 & 3 capacity 12 plates high x		
12 = 144 x 2 shelves	=	288 dishes
Total shelf capacity	=	348 dishes

NOTE: If the area under the shelf assembly is to be used for additional plate storage, then add 84 more plates.

Suspended shelves for #830-SERIES glove boxes.

Catalog #800-SHELF-830 SUSPENDED SHELF.


The total number of shelves for the **#830-SERIES** glove boxes will vary based upon the overall size of the chamber selected. Once you have selected the #830 chamber size you want, contact us about the corresponding shelf size to fit.

The photo on left is a typical set up for the **#830-ABC** Compact Glove box. Each shelf has a capacity of fifty (50) Petri Dishes (plates), for a total 3 shelf capacity of 150. If plates are to be stored below the shelves, there is room for an additional forty (40) plates, making a total of 190 plates.

Accessories for all Chambers and Glove Boxes

ENTRY PORT VARIATION

PLAS ■ LABS offers numerous variations for access to the internal environment of a glove box or isolator.



1. #800-IRIS/S

REPLACEMENT IRIS SLEEVES and CUFFS. Refer to **#800-IRIS** for description.

2. #800-IRIS

IRIS DIAPHRAGMS (UNMOUNTED). A flexible cuff is rotated to tighten around the operators arms. Opening diameter is 6" (152 mm) I.D. Sold in pairs only.

3. #800-IRIS/CONV

IRIS DIAPHRAGMS CONVERSION SET. Iris diaphragms are mounted on formed plastic caps that fit over our standard glove port rings. Twist into place. Sold in pairs only.

4. #800-BH/CONV

BARE HAND ENTRY PORT CONVERSION. These ports are 8"Ø (203 mm) in diameter and consist of five (5) pleats of flexible Hypalon™. Sold in pairs only. Ten (10"Ø) diameter (524 mm) ports are also available.

Did you Know?
These Bare Hand entries can be mounted on any glove box or isolator.

NOTE: #800-IRIS/MOUNT

FACTORY MOUNTED IRIS DIAPHRAGMS. Refer to **#800-IRIS** for description.

GLOVE PORT PLUGS



#800-PLUG & #815-PLUG (La Petite Unit)

Plug Applications

#800-PLUG & #815-PLUG Glove Port Plugs are used for sealing ports during:

- Glove change over.
- Bare-hand entry change over.
- Iris diaphragm change over.
- When ports are not needed.



Insert plug into port opening and rotate red knob (tighten) until port is sealed.

Accessories for all Chambers and Glove Boxes

ATMOSPHERE PURIFICATION & STERILIZATION ITEMS

When your research requires inner atmosphere purification and/or sterilization, we have the answer for your needs. Particle filtration, harmful vapors, or just static removal devices are readily available for your glove box or isolators.

Our sterilants and cleaners are all F.D.A. approved.

There are several commercially available sterilizing agents and technologies that can be used on our cast acrylic isolators and glove boxes. These include:

- CLIDOX-S™ Sterilant by Pharmalabs™.
- Chlorine Dioxide Gas Tablets (MB-10 Tablets by QUIP Labs™).
- Formaldehyde.
- Diluted isopropyl alcohol. 70% (v/v) in water (Bactericide).
- Diluted ethyl alcohol (Virucide).
- Chlorine (sodium hypochlorite) 5 g/Liter

NOTES about our new Carbon filter types and set up.



Typical Carbon Filter for glove boxes and isolators. Refer to notes this page.



Typical Carbon Filter for glove boxes and isolators. Refer to notes this page.

Planning on using "VHP" sterilization equipment?



#818-GB/SPCL for use with a VHP sterilizer

NOTE: "VHP" sterilization requires a special plastic, such as clear PVC or PETG. "Banjo™" style purge fittings are included with the glove box.



Close up of "BANJO™" style purge fitting Catalog #800-BANJO. Refer to page #46.

Carbon Filter Specifications and Types.

- **#800-HEPA/CARB.** Combination HEPA and Non-Impregnated Carbon absorbs most aliphatic and aromatic hydrocarbons, solvents, organic vapors, aldehydes, ketones, alcohols, organic acids, esters, sulfur compounds, and odors.
- **#800-CARB (A).** Impregnated Carbon filter for alkaline type fumes.
- **#800-CARB (B).** Impregnated Carbon filter for ammonia or amines.
- **#800-CARB (C).** Non-Impregnated Carbon filter for absorption of aromatic hydrocarbons, organic vapors, ketones, alcohols, organics acids, esters, Nitrogen compounds, and odors.

Accessories for all Chambers and Glove Boxes

ATMOSPHERE PURIFICATION (CONTINUED)




#800-HEATER & #800-PC
Specify Export or Domestic

#800-HEATER (1.), Catalyst Heater Unit

Maintains correct incubation temperature levels and helps reduce trace amounts of Oxygen in the glove box. Includes **thermocouple (2.)** for temperature sensing. Place thermocouple at your desired height inside chamber.

#800-PC (3.) Palladium Canister filled with 1/2 Lb. of Palladium (Pd). Canisters sold separately.

The heater draws the internal atmosphere off the floor and pushes it up across heating elements and Palladium pellets. Any trace amounts of Oxygen are then reduced into water (H_2O) vapor. The Drying Train **#800-DT** (see below) then removes the vapor as needed.



#800-ONEG Oxygen Removal System
#800-PUMP Pump not included

#800-ONEG includes O_2 indicator (1.) Oxygen Trap (2.) hoses and plastic holder.

#800-ONEG reduces O_2 level down to 10-12 ppb. Must order as a set only.

NOTE: Glove box must be made anaerobic before attaching this system.



#800-DT Drying Train with Pump
Specify Export or Domestic

#800-DT removes unwanted moisture (H_2O) from the internal atmosphere of the glove box. Includes pump, hoses, and three (3) canisters filled with Molecular Sieve™ (**#800-DESI**).

#800-DESI includes three (3) canisters filled with Molecular Sieve™ and hoses with quick disconnects. Pump not included.



#800-PUMP (Diaphragm type)
Specify Export or Domestic.

#800-PUMP, Diaphragm pump mounted on a .375" plastic base with suction cup feet for stability.

Includes an in-line on/off switch and five (5) feet of flexible hose. No quick disconnects.

Electrical requirements for all products this page:

North American: 115-120 Volts 60 Hz. 5 Amps.

International (Export): 220-240 Volts 50 Hz. 2.5 Amps.

Additional Items for Filtration.

- #CF-50/9** Carbon filter pad, 9"Ø (228 mm). "Spun Media" color black. Twelve (12) per package.
- #GF-50-B9** Glass media Filter Pad 9"Ø (228 mm) Color pink, (5 μ) microns. Twelve (12) per package.
- #MF-50/9** Spunbound Polypropylene pre-filter, 9"Ø (228 mm), color white, 0.4 μ non-hygroscopic. Twelve (12) per package.

Accessories for all Chambers and Glove Boxes

STAINLESS STEEL SUPPORT CARTS (TROLLEYS)



#800- Cart III

These carts are very useful for increasing portability and laboratory space.

All carts have anti-scuff casters including brakes and locks. Carts also feature a cushioning system that minimizes vibration, shocks, and wear and tear on chambers.

Most carts feature horizontal gas tank support racks which provide greater tank handling safety and flexibility.

Also available in a “knocked down” (unassembled) version. See information below.



#800-CART GB
with gas tank support racks

PRODUCT NUMBER	CART DESCRIPTION	SIZES (O.D.)
CART-GB	For #GB, #NB, #AC Units	51”w x 24”d x 35” h
CART-GBB	For #GBB, #NBB, # ACB Units	55”w x 34”d x 35” h
CART-ABC	For #ABC Glove Box	43”w x 24”d x 34” h
CART-ABB	For #ABB Glove Box	49”w x 30”d x 34” h
CART-ABD	For #ABD Glove Box	63”w x 31”d x 34” h
CART-ABE	For #ABE Glove Box	76”w x 32”d x 34” h
CART-CLC	For #870 sized units	55”w x 34”d x 31” h

Custom Sizes and Designs available



Kit includes: Main base frame, casters, uprights, and top.



Insert casters and tighten.



Force fit the four legs (uprights).



Place top on four legs & press into place.



#800- Cart KD

SAVES SHIPPING COSTS!

STAINLESS STEEL “KNOCKED DOWN” CART (TROLLEY)

The Knocked Down Cart (Trolley) size is 40” (1016 mm) w x 24” (610 mm) d x 34” (864 mm) h.


Accessories for all Chambers and Glove Boxes

MICROSCOPE & INSPECTION LENS COMPONENTS

Visual inspection of your research components can be made without removing them from the isolator or glove box.

Our new inspection lens and microscope port accessories can be factory mounted on the viewing surface of any unit we manufacture.

Refer to the photo on the right.



#800-MICRO and #800-BH/CONV

#800-MICRO MICROSCOPE VIEWING PORT is 6.25" Ø (159 mm) and includes the port, Hypalon™ cover, and stainless steel worm clamp. The Hypalon™ cover allows the microscope oculars to protrude from the chamber to permit exterior viewing of your research materials. Specify single or double oculars.

#800-HYCOV REPLACEMENT COVER for the **#800-MICRO**. Sold in pairs only. One with ocular hole(s) and one without for use when the microscope is not needed. Specify your preference when ordering.

#800-BH/CONV BARE HAND CONVERSION KIT.
Refer to page #39 for complete details.

NEW ITEM!



#800-I/LENS (Inspection Lens)

The new inspection lens is rated having a focal length tolerance of $\pm 1\%$ at 11.8" inches (300 mm). The diameter of the lens is 6" (152 mm). This is a "fused silica" glass lens of extremely high optical quality. Factory mounted only.



#800-STAND

The **#800-STAND** has an adjustment span of 9" (228 mm) and is 16.5" (420 mm) high. The base has non-slip foot pads and is 8" x 10" (203 x 254 mm).

Accessories for all Chambers and Glove Boxes

GLOVE BOX & ISOLATOR GLOVES

Material selections.

PLAS ■ LABS offers one piece gauntlet style gloves in four (4) different materials. They are:

1. Butyl™ **Catalog #800-GB-Series** (Not shown).
2. Hypalon™ **Catalog #800-GH-Series.**
3. Latex **Catalog #800-GLLP-Series.**
4. Neoprene™ **Catalog #800-GN-Series.**

Refer to photo on right.

We also offer glove box sleeves with replaceable hands. The four (4) different hand materials are:

1. Butyl™ **Catalog #800-GA-BUT-Series.**
2. Hypalon™ **Catalog #800-GA-HYP-Series.**
3. Neoprene™ **Catalog #800-GA-NEO-Series.**
4. Nitrile™ **Catalog #800-GA-NIT-Series.**

Refer to photo below.



One piece gauntlet style gloves



Replaceable hands for gauntlet sleeves

NOTES: The #800-GH Hypalon™ glove is the standard glove used on most Plas ■ Labs products.

The #800-GLLP (latex glove) is the primary glove used on the La Petite glovebox (Catalog #815-PG).

Some port diameters are available in three (3) sizes ranging from 7" Ø (178 mm) to 10" Ø (254 mm).


Hand sizes are also available varying in size from 7" Ø (177 mm) to 10" Ø (254 mm).

Refer also to the gloves/hands data chart for additional details. Refer to page #45.

Glove box sleeve with replaceable hand set up. Components are:

1. Typical sleeve with mounting cuff. **Catalog #800-GH/S.**
2. "O" rings for sealing glove to sleeve.
3. Hand (glove) replacement. **Catalog #800-GA/NIT.**

Note: Nitrile hand component illustrated (Both photos). Other hand materials available.



Replaceable hand set up with gauntlet sleeve



Assembled hand/sleeve set up

Accessories for all Chambers and Glove Boxes

GLOVE BOX & ISOLATOR GLOVE SPECIFICATIONS

We offer 22 variations of gauntlet gloves and replacement hands for glove boxes and isolators. All of the gauntlet gloves (sleeves) can also be mounted onto our conversion kits. Conversion kits (formed plastic caps) can be mounted onto any of our glove boxes or isolators. [Refer also to page #39.](#)

Specification Chart

Product Number	Description	Material	Port Diameter	Hand Size	Thickness
800-GH/S	Sleeve Only	Hypalon	8"	n/a	.015"
800-GOS	Glove over sleeve connector				
Optional hand replacement for use with 800-GH/S, sleeve only. 800-GOS, glove over sleeve connector, is needed for use with the hand replacements and sleeves					
800-GA/BUT/7	Hand Replacement	Butyl	n/a	7"	.015"
800-GA/BUT/8	Hand Replacement	Butyl	n/a	8"	.015"
800-GA/BUT/9	Hand Replacement	Butyl	n/a	9"	.015"
800-GA/BUT/10	Hand Replacement	Butyl	n/a	10"	.015"
800-GA/HYP/7	Hand Replacement	Hypalon	n/a	7"	.015"
800-GA/HYP/8	Hand Replacement	Hypalon	n/a	8"	.015"
800-GA/HYP/9	Hand Replacement	Hypalon	n/a	9"	.015"
800-GA/HYP/10	Hand Replacement	Hypalon	n/a	10"	.015"
800-GA/NEO/7	Hand Replacement	Neoprene	n/a	7"	.015"
800-GA/NEO/8	Hand Replacement	Neoprene	n/a	8"	.015"
800-GA/NEO/9	Hand Replacement	Neoprene	n/a	9"	.015"
800-GA/NEO/10	Hand Replacement	Neoprene	n/a	10"	.015"
800-GA/NIT/7	Hand Replacement	Nitrile	n/a	7"	.015"
800-GA/NIT/8	Hand Replacement	Nitrile	n/a	8"	.015"
800-GA/NIT/9	Hand Replacement	Nitrile	n/a	9"	.015"
800-GA/NIT/10	Hand Replacement	Nitrile	n/a	10"	.015"
One Piece Gloves					
800-GB	One-Piece Glove	Butyl	8"	9"	.015"
800-GB/7	One-Piece Glove	Butyl	8"	7"	.015"
800-GB/8	One-Piece Glove	Butyl	8"	8"	.015"
800-GB/9/6	One-Piece Glove	Butyl	7"	9"	.015"
800-GB/9/10	One-Piece Glove	Butyl	10"	9"	.015"
800-GH	One-Piece Glove	Hypalon	8"	9"	.015"
800-GH/7	Two-Piece Glove	Hypalon	8"	7"	.015"
800-GH/8	One-Piece Glove	Hypalon	8"	8"	.015"
800-GH/10	One-Piece Glove	Hypalon	8"	10"	.015"
800-GH/8/10	One-Piece Glove	Hypalon	10"	8"	.015"
800-GH/9/10	One-Piece Glove	Hypalon	10"	9"	.015"
800-GH/9/6	One-Piece Glove	Hypalon	7"	9"	.015"
815-PGB, La Petite Glove Box Gloves					
800-GLLP	One-Piece Glove	Latex	7"	9"	.015"
800-GLLP/7	One-Piece Glove	Latex	7"	7"	.015"
800-GLLP/10	One-Piece Glove	Latex	7"	10"	.015"
800-GN	One-Piece Glove	Neoprene	8"	9"	.015"
800-GN/8	One-Piece Glove	Neoprene	8"	8"	.015"
800-GN/10	One-Piece Glove	Neoprene	8"	10"	.015"
800-GN/11	One-Piece Glove	Neoprene	8"	11"	.015"
800-GN/9/6	One-Piece Glove	Neoprene	7"	9"	.015"

Accessories for all Chambers and Glove Boxes

VALVE AND GAUGE VARIATIONS

Various valves, gauges, and meters can be added to any glove box, isolator, or desiccator. These are normally mounted in standard locations, but we can easily change set ups depending upon your own preferences.

1. **#800-IN/OUT. COMPLETE VALVE PACKAGE** (Unmounted). Includes ground key cock valve, with flanges, and inner/outer serrated nipples.
2. **#800-GOFAU. GOOSE NECK WATER FAUCET.** Factory mounted internal faucet has a swivel goose neck for convenience. Useful for controlling water supply into the chamber. Includes serrated nipple and flange for easy hose connection.
3. **#800-FLOW. FLOW METER KIT.** Helps regulate gas flow at lower incoming pressures. Kit includes flowmeter, mounting bracket, hoses and ring clips. Range "0 -10 LPM." Refer to photo below also (#3).
4. **#800-HYGRO. HYGROMETER** (Large). Indicates inner level of chamber Relative Humidity (Rh).




5. **#800-FLOW/D. SMALL FLOW METER** used for smaller chambers and desiccators. Helps regulate gas flow at lower incoming pressures. Kit includes flowmeter, mounting bracket, hoses and ring clips. Range "5-50 SCFH."
6. **#800-BANJO. BANJO™ TYPE FITTINGS** for use with VHP sterilization equipment and procedures. O.D. is 1.375" (35 mm). I.D. is .875" (23 mm).
7. **#800-NLP. SERRATED NIPPLE ONLY.** Normally attached to the **#800-VALVE** on the inside of the chamber.
8. **#800-VALVE. GROUND KEY COCK VALVE** (Only). Valve with flange can also be mounted on a 90° elbow if desired. Specify N₂ or VAC (Vacuum) valve label.
9. **#800-PRV. PRESSURE RELIEF VALVE.** Minimizes the chance for chamber implosion or explosion because of incorrect pressure level.

#800-MAGNA. Magnehelic Gauge Package.

These Magnehelic gauges indicate the internal pressure levels of the HEPA filter system and the internal pressure of the glove box.

- A. This **RED** pointer indicates when it is time to change the HEPA filters.
- B. When the **RED** pointer lines up with this "set pointer," it is time to change the HEPA filters.
- C. This **RED** pointer will change because of HEPA filter variations in resistance.
- D. This is the main electrical power (On/Off) switch.



#800-MAGNA Magnehelic Gauge Package (H₂O)

Accessories for all Chambers and Glove Boxes

MISCELLANEOUS ACCESSORIES

Numerous systems are available for gaining access into any glove box or isolator. Some of these can be set up so as not to disturb the interior atmosphere of the chamber.



#800-BLANK/P BLANKING PLATE PORT

#800-BLANK/P BLANKING PLATE PORT. Useful for attaching special cables and plugs. This kit includes a blank aluminum plate, fastening stainless steel nuts, bolts, and sealant. The blank plate itself is ready for machining as necessary by the end user. The hole cut out is 6" x 3" (152 x 76 mm) and the blank plate is 4.5" x 7.5" (190 x 114 mm).

NOTE: Other sizes can be requested (Special quote).




#800-CMPTR COMPUTER INTERFACE KIT

Permits interface connection between internal devices and your external computer. Kit includes stress relief plates, shielded computer cable (8 feet), stainless steel fasteners, and RS-232 connectors.


NOTE: Other connectors can be supplied.

*THE NEW "INTERLOCK SYSTEM"
for Transfer Chamber doors.*



#800-DOOR/ILOCK AUTOMATIC DOOR LOCK SYSTEM

Exterior door open and inner door locked



#800-DOOR/ILOCK AUTOMATIC DOOR LOCK SYSTEM

Exterior door open and inner door locked

This automatic door locking system prevents both transfer chamber doors from being opened at the same time. The interlock system does not work unless one door is opened. If one door is opened, the other door is locked in five (5) seconds. The locking mechanism is rated at 120 Lbs of force. This is for operator safety. It operates on a 12 volt system. Power supply is included.

Electrical requirements for this system:

North American: 115-120 Volts 60 Hz. 1 amp.

International (Export): 220-240 Volts 50 Hz. 1 amp.

Accessories for all Chambers and Glove Boxes

ADDITIONAL MISCELLANEOUS ACCESSORIES


"No project is too small !"

PLAS ■ LABS offers a variety of unique miscellaneous accessories for our glove boxes, isolators, and chambers. All accessories are easily adaptable to standard or customized research equipment.

If you are doing research requiring unique enclosures, please contact us. We might be able to assist you in your research efforts.

"New Dip Tank"

#800-TANK, GERMICIDAL DIP TANK. Useful for quick entry into and out of the chamber. A "deep draw" bottom section of the chamber is required for this accessory. Package includes inner and outer tank with stainless steel exterior lid. Fill the tank with germicidal solution to sanitize items being introduced into the chamber.



#800-TANK, GERMICIDAL (DIP) DUNK TANK


#800-GTREG, Gas Tank Regulator. Regulator features two stainless steel diaphragms that will not contaminate the gas stream and is resistant to in-board diffusion. The design also features low dead volume (internal volume), which minimizes particle generation. Kit includes 5 feet of hose, (.250 I.D. x .500 O.D.) and attachment clamp.



#800-GTREG

Very useful for scaling down high pressure levels prior to reaching flow meters. Domestic (US) thread styles only.

#800-T, "T" SHAPED TRANSFER CHAMBER. This is used for connecting two (2) glove boxes together. Included in the kit are a leveling support stand, one (1) extra vacuum gauge, two (2) extra gas key cock valves, two (2) leveling trays for liquids, and connecting tape. The center section has its own access door with stainless steel locking clamps.



#800-T, "T" SHAPED TRANSFER CHAMBER



#800-02, OXYGEN INDICATOR-ICS

#800-02, OXYGEN INDICATOR. A small hand held device to monitor Oxygen levels inside the chamber. It is battery operated and has an effective reading range of 0-100%.

PLAS ■ LABS™



*Typical training session
in the Plas ■ Labs "Wet Lab"*

Good training is essential to be prepared for future business opportunities. Our training facility can accommodate up to six (6) people at a time. A typical training session takes about two (2) days, and Plas ■ Labs pays all hotel, and meal (entertainment) expenses during that time.

Visit us and gain the knowledge and confidence to prepare you for success when selling our products. After training, you will quickly start counting your rewards.

PLAS ■ LABS™

401 E. NORTH ST.
LANSING, MICHIGAN
48906 U.S.A.

Innovated Quality Since 1967

www.plas-labs.com ■ E-mail: sales@plas-labs.com ■ 1.517.372.7177 ■ 1.800.866.7527 ■ Fax: 1.517.372.2857

PLAS LABS™

Clearly Your Best Choice™



"We Produce Equipment That Helps You Produce Results"

- Anaerobic Chambers
- Cadaver / Prosection Trays
- Tissue Culture Hoods
- Desiccators
- Stress Simulators
- Earth Science Equipment
- Sedimentation Chambers
- Dust Free Hoods
- Demonstration Hoods
- Controlled Atmosphere Chambers
- Guinea Pig Housing Units
- Isolation Chambers
- Custom Fabrication
- Intensive Care Units
- Oxygen Therapy Units
- Rabbit Restainers
- Rodent Restainers
- Primate Chairs
- Small Animal Surgery Boards
- Gnotobiotic Isolators
- Glove Boxes
- Metabolism Chambers
- Euthanasia Chambers
- Balance Chambers
- Animal "Ward" Housing Units
- Canine Kennel Runs

DISTRIBUTED BY:

PLAS ■ LABS™

401 E. NORTH ST.
LANSING, MICHIGAN
48906 U.S.A.

Innovated Quality Since 1967

www.plas-labs.com ■ E-mail: sales@plas-labs.com ■ 1.517.372.7177 ■ 1.800.866.7527 ■ Fax: 1.517.372.2857