




The name you can trust

Innovation | Technology | Precision



Humidity (Stability) Chambers


Cooling Chambers & Deep Freezers

Incubators

Ovens & Vacuum Ovens


Web Based Software

Web Based Data Logging System



Newtronic Lifecare Equipment Pvt. Ltd.
www.newtronic.in

Introducing Remote View & Alarm System(RVAS)



RVAS

Remote View & Alarm System


Take Remote Monitoring to an in-depth level of information-sharing with NEWTRONIC's Remote View and Alarm System (RVAS).

RVAS Features:

1. View & Acknowledge real-time T/RH readings of EACH sensor
2. View status of each I/O Components (humidifier, heater, compressor, etc)
3. View status of Utilities (Mains Supply, Door, Water Level)
4. Generate SMS for critical alarms. SMS alerts to specific mobile nos. based on equipment ownership
5. Get real-time video feeds from IP-cameras



Remote View & Alarm System



Other Features:

- 7" colour touch screen display
- Connect upto 20 equipments on TCP/IP with on-site programming
- Eliminate the need for cabling
- Audio hooter & visual display to notify alarms & deviations
- Department-wise segregation of users & chambers
- System Diagnosis
- IQ, OQ protocols



Real Time Status



SMS Alerts



Alerts 24/7



Alert Groups



Multiple Interfaces



IP Camera

Touching tomorrow....



Constant focus on technology & innovation

Leading in diversity

Strong culture of performance

Believing in trust & teamwork

Ambitious long term goals

Sustainable Strategy

NL/CC/PH/0514/1.2

Newtronic Service centers:

Domestic:

Ahmedabad	Hyderabad / Vizag
Baddi	Indore
Bangalore	Kolkata / Sikkim
Chennai	Mumbai
Delhi	Pune
Goa	

International:

Australia	Mexico
Bangladesh	Nepal
China	New Zealand
Dominican Republic	Phillipines
Egypt	Singapore
Indonesia	Thailand
Jordan	Turkey
Malaysia	

service@newtronic.in

service.exports@newtronic.in

Innovation | Technology | Precision



The name you can trust

NEWTRONIC LIFECARE EQUIPMENT PVT. LTD. ISO 9001:2008 Certified Company

108-ABCD, Kandivali Co-op. Ind. Estate Ltd.
Charkop, Kandivli (West),
Mumbai - 400067, India
T: +91 -22- 2867 9326 |
F: +91 -22- 2867 6059

Email:

India : sales@newtronic.in

International : exports@newtronic.in

Visit us on: www.newtronic.in