

EN

Disposable or reusable

/// Autoclavable disposable tubes for the MultiDrive control

Avoiding cross-contamination effectively: The BT 250.5 and MT 150.5 disposable tubes for the MultiDrive control are designed to be autoclavable and reusable, too.

FOR WET SAMPLES

Blending tube

BT 250.5 is a disposable tube with a cross-shaped stainless steel beater and flow breakers.

It is ideally suited for wet blending.



FOR DRY SAMPLES

Mixing tube

MT 150.5 is a disposable tube with a typical stainless steel beater.

It is ideally suited for grinding dry samples.



/// Good to know

- › Both tubes are autoclaveable and thus prevent cross-contamination effectively even when used multiple times.
- › Spare seals extend the service life of both tubes.
- › An integrated RFID chip ensures the perfect lifespan. You will get informed when the seal or the tube need to be changed.
- › MultiDrive control recognizes automatically, if the tube is covered by the mandatory protective housing – otherwise it won't start.

TECHNICAL DATA

	BT 250.5	MT 150.5
Quantity	5 tubes and 5 spare seals per package; reactivation codes included	
Useable volume max.	250 ml	150 ml
Feed grain size / Feed hardness max.	10 mm / 5 Mohs	
Material in contact with medium	Stainless steel, PCTG	
Autoclavable	Yes	
FDA compliance	Yes	

Only usable in combination with TC 1 protective housing. TC 1 is available as an accessory.

